



Gelest

Silane Coupling Agents: *Connecting Across Boundaries*

Enhance Adhesion

Increase Mechanical Properties

Improve Dispersion

Provide Crosslinking

Immobilize Catalysts

Bind Biomaterials

Version 2.0:

*New Coupling Agents
for Metal Substrates !*

*New Coupling Agents for Vapor Phase
Deposition !*

*New Coupling Agents
for Proteins !*



Gelest, Inc.

Telephone: General 215-547-1015
Order Entry 888-734-8344
FAX: 215-547-2484
Internet: www.gelest.com
Correspondence:
11 East Steel Road
Morrisville, PA 19067, USA

In Europe: ABCR GmbH & Co. KG
Im Schleht
D-76187 Karlsruhe
Germany
Tel: +49 - 721 - 950610
Fax: +49 - 721 - 9506180
e-mail: info@abcr.de
on-line catalog: www.abcr.de

In Japan: AZmax Co. Ltd. Tokyo Office
Matsuda Yaesudori Bld F8
1-10-7 Hatchoubori, Chuo-Ku
Tokyo 104-0032
Tel: 81-3-5543-1630
Fax: 81-3-5543-0312
email: sales@azmax.co.jp
on-line catalog: www.azmax.co.jp

In South-East Asia:

Altus Technologies Pte Ltd
196 Pandan Loop #06-09
Pantech Industrial Complex Singapore 128384
Tel: (65) 6779 7666 Fax: (65) 6779-7555
www.altus.com.sg

For further information consult our website at: www.gelest.com

Sales of all products listed are subject to the published terms and conditions of Gelest, Inc.

Gelest **Silane Coupling Agents**

Connecting Across Boundaries



TABLE OF CONTENTS

What is a Silane Coupling Agent?	2
How Does a Silane Coupling Agent Work?	3
Selecting a Silane Coupling Agent - Inorganic Substrate Perspective	4
Selecting a Silane Coupling Agent - Polymer Applications	5
Selecting a Silane Coupling Agent - Interphase Considerations	9
Special Topics:	
Dipodal Silanes	11
Linker Length	12
Cyclic Azasilanes	13
Thermal Stability of Silanes	14
Aqueous Systems & Water-Borne Silanes	15
Masked Silanes - Latent Functionality	16
Coupling Agents for Metal Substrates	17
Difficult Substrates	18
Applying a Silane Coupling Agent	19
Silane Coupling Agents for Polymers - Selection Chart	21
Silane Coupling Agents for Biomaterials - Selection Chart	24
Silane Coupling Agents - Properties	25
Organosilane-Modified Silica Nanoparticles	54
Further Information - Other Resources	55
Index	56

Silane Coupling Agents: Connecting Across Boundaries v2.0
by Barry Arkles




What is a Silane Coupling Agent?

Silane coupling agents have the ability to form a durable bond between organic and inorganic materials. Encounters between dissimilar materials often involve at least one member that's siliceous or has surface chemistry with siliceous properties; silicates, aluminates, borates, etc., are the principal components of the earth's crust. Interfaces involving such materials have become a dynamic area of chemistry in which surfaces have been modified in order to generate desired heterogeneous environments or to incorporate the bulk properties of different phases into a uniform composite structure.



The general formula for a silane coupling agent typically shows the two classes of functionality. X is a hydrolyzable group typically alkoxy, acyloxy, halogen or amine. Following hydrolysis, a reactive silanol group is formed, which can condense with other silanol groups, for example, those on the surface of siliceous fillers, to form siloxane linkages. Stable condensation products are also formed with other oxides such as those of aluminum, zirconium, tin, titanium, and nickel. Less stable bonds are formed with oxides of boron, iron, and carbon. Alkali metal oxides and carbonates do not form stable bonds with Si-O-. The R group is a nonhydrolyzable organic radical that may possess a functionality that imparts desired characteristics.

The final result of reacting an organosilane with a substrate ranges from altering the wetting or adhesion characteristics of the substrate, utilizing the substrate to catalyze chemical transformations at the heterogeneous interface, ordering the interfacial region, and modifying its partition characteristics. Significantly, it includes the ability to effect a covalent bond between organic and inorganic materials.



How does a Silane Coupling Agent Work?

Most of the widely used organosilanes have one organic substituent and three hydrolyzable substituents. In the vast majority of surface treatment applications, the alkoxy groups of the tri-alkoxysilanes are hydrolyzed to form silanol-containing species. Reaction of these silanes involves four steps. Initially, hydrolysis of the three labile groups occurs. Condensation to oligomers follows. The oligomers then hydrogen bond with OH groups of the substrate. Finally during drying or curing, a covalent linkage is formed with the substrate with concomitant loss of water. Although described sequentially, these reactions can occur simultaneously after the initial hydrolysis step. At the interface, there is usually only one bond from each silicon of the organosilane to the substrate surface. The two remaining silanol groups are present either in condensed or free form. The R group remains available for covalent reaction or physical interaction with other phases.

Silanes can modify surfaces under anhydrous conditions consistent with monolayer and vapor phase deposition requirements. Extended reaction times (4-12 hours) at elevated temperatures (50°-120°C) are typical. Of the alkoxy silanes, only methoxysilanes are effective without catalysis. The most effective silanes for vapor phase deposition are cyclic azasilanes.

Hydrolysis Considerations


Water for hydrolysis may come from several sources. It may be added, it may be present on the substrate surface, or it may come from the atmosphere.

The degree of polymerization of the silanes is determined by the amount of water available and the organic substituent. If the silane is added to water and has low solubility, a high degree of polymerization is favored. Multiple organic substitution, particularly if phenyl or tertiary butyl groups are involved, favors formation of stable monomeric silanols.


The thickness of a polysiloxane layer is also determined by the concentration of the siloxane solution. Although a monolayer is generally desired, multilayer adsorption results from solutions customarily used. It has been calculated that deposition from a 0.25% silane solution onto glass could result in three to eight molecular layers. These multilayers could be either interconnected through a loose network structure, or intermixed, or both, and are, in fact, formed by most deposition techniques. The orientation of functional groups is generally horizontal, but not necessarily planar, on the surface of the substrate.

The formation of covalent bonds to the surface proceeds with a certain amount of reversibility. As water is removed generally by heating to 120°C for 30 to 90 minutes or evacuation for 2 to 6 hours, bonds may form, break, and reform to relieve internal stress. The same mechanism can permit a positional displacement of interface components.

Hydrolytic Deposition of Silanes



Anhydrous Deposition of Silanes



Selecting A Silane Coupling Agent - Inorganic Substrate Perspective

Factors influencing silane coupling agent selection include:


- Concentration of surface hydroxyl groups*
- Type of surface hydroxyl groups*
- Hydrolytic Stability of the bond formed*
- Physical dimensions of the substrate or substrate features*

Coupling is maximized when silanes react with the substrate surface and present the maximum number of sites with reactivity specific and accessible to the matrix phase. An additional consideration is the physical and chemical properties of the interphase region. The interphase can promote or detract from total system properties depending on its physical properties such as modulus or chemical properties such as water/hydroxyl content.

Hydroxyl-containing substrates vary widely in concentration and type of hydroxyl groups present. Freshly fused substrates stored under neutral conditions have a minimum number of hydroxyls. Hydrolytically derived oxides aged in moist air have significant amounts of physically adsorbed water which can interfere with coupling. Hydrogen bonded vicinal silanols react more readily with silane coupling agents, while isolated or free hydroxyls react reluctantly.

Silane coupling agents with three alkoxy groups are the usual starting point for substrate modification. These materials tend to deposit as polymeric films, effecting total coverage and maximizing the presentation of organic functionality. They are the primary materials utilized in composites, adhesives, sealants, and coatings. Limitations intrinsic in the utilization of a polylayer deposition are significant for nano-particles or nano-composites where the interphase dimensions generated by polylayer deposition may approach those of the substrate. Residual (non-condensed) hydroxyl groups from alkoxy-silanes can also interfere in activity. Monoalkoxy-silanes provide a frequently used alternative for nano-featured substrates since deposition is limited to a monolayer.

If the hydrolytic stability of the oxane bond between the silane and the substrate is poor or the application is an aggressive aqueous environment, dipodal silanes often exhibit substantial performance improvements. These materials form tighter networks and may offer up to 10⁵x greater hydrolysis resistance making them particularly appropriate for primer applications.



Amino-silanes couple fiberglass to phenolic or urea-formaldehyde resins

Silane Effectiveness on Inorganics

	SUBSTRATES
EXCELLENT ↑ GOOD ↑ SLIGHT ↑ POOR	Silica
	Quartz
	Glass
	Aluminum (AlO(OH))
	Alumino-silicates (e.g. clays)
	Silicon
	Copper
	Tin (SnO)
	Talc
	Inorganic Oxides (e.g. Fe ₂ O ₃ , TiO ₂ , Cr ₂ O ₃)
POOR	Steel, Iron
	Asbestos
	Nickel
	Zinc
	Lead
	Marble, Chalk (CaCO ₃)
	Gypsum (CaSO ₄)
	Barytes (BaSO ₄)
Graphite	
Carbon Black	

Estimates for Silane Loading on Siliceous Fillers

Average Particle Size	Amount of Silane (minimum of monolayer coverage)
<1 micron	1.5%
1-10 microns	1.0%
10-20 microns	0.75%
>100 microns	0.1% or less


Selecting A Silane Coupling Agent - Polymer Applications

Coupling agents find their largest application in the area of polymers. Since any silane that enhances the adhesion of a polymer is often termed a coupling agent, regardless of whether or not a covalent bond is formed, the definition becomes vague. In this discussion, the parochial outlook will be adopted, and only silanes that form covalent bonds directly to the polymer will be considered. The covalent bond may be formed by reaction with the finished polymer or copolymerized with the monomer. Thermoplastic bonding is achieved through both routes, although principally the former. Thermosets are almost entirely limited to the latter. The mechanism and performance of silane coupling agents is best discussed with reference to specific systems. The most important substrate is E-type fiberglass, which has 6-15 silanol groups per μm^2 .


Thermosets

Acrylates, methacrylates and Unsaturated Polyesters, owing to their facility for undergoing free-radical polymerization, can be modified by copolymerization with silanes that have unsaturated organic substitution. The usual coupling agents for thermoset polyesters undergo radical copolymerization in such systems. These resins, usually of loosely defined structure, often have had their viscosity reduced by addition of a second monomer, typically styrene. In general, better reinforcement is obtained when the silane monomer matches the reactivity of the styrene rather than the maleate portion of the polyester.


Methacrylyl and styryl functional silanes undergo addition much more readily than vinylsilanes. A direct approach to selecting the optimum silane uses the e and Q parameters of the Alfrey-Price treatment of polymerization. Here e indicates the polarity of the monomer radical that forms at the end of a growing chain, while Q represents the resonance stabilization of a radical by adjacent groups. Optimum random copolymerization is obtained from monomers with similar orders of reactivity. Vinyl functional silanes mismatch the reactionary parameters of most unsaturated polyesters. However, they can be used in direct high pressure polymerization with olefins such as ethylene, propylene and dienes.



Acrylate Coupling Reaction



Unsaturated Polyester (Styrene) Coupling Reaction



Polyethylene Graft Coupling Reaction

Urethanes

Thermoset urethane can be effectively coupled with two types of silanes. The first type, including isocyanate functional silanes, may be used to treat the filler directly or integrally blended with the diisocyanate (TDI, MDI, etc.) prior to cure. Amine and alkanolamine functional silanes, on the other hand, are blended with the polyol rather than the diisocyanate. Isocyanate functional silanes couple with the polyol. Alkanolamine functional silanes react with the isocyanate to form urethane linkages, while amine silanes react with the isocyanates to yield urea linkages. A typical application for coupled urethane system is improving bond strength with sand in abrasion-resistant, sand-filled flooring resins.

Moisture-Cureable Urethanes


Aminosilanes have the general ability to convert isocyanate functional urethane prepolymers to systems that crosslink in the presence of water and a tin catalyst. The preferred aminosilanes are secondary containing methyl, ethyl or butyl substitutions on nitrogen.

Epoxies


Epoxycyclohexyl and glycidoxy functional silanes are used to pretreat the filler or blended with the glycidylbisphenol-A ether. Amine functional silanes can likewise be used to pretreat the filler or blended with the hardener portion. Treatment of fillers in epoxy adhesives improves their dispersibility and increases the mechanical properties of the cured resin. A large application area is glass cloth-reinforced epoxy laminates and preregs in aerospace and electrical printed circuit board applications.

Phenolics


Phenolic resins are divided into base catalyzed single-step resins called resols or better known acid catalyzed two-step systems called novolaks. Although foundry and molds are formulated with resols such as aminopropylmethyl-dialkoxysilanes, the commercial utilization of silanes in phenolic resins is largely limited to novolak/glass fabric laminates and molding compounds. The phenolic hydroxyl group of the resins readily react with the oxirane ring of epoxy silanes to form phenyl ether linkages. When phenolic resins are compounded with rubbers, as in the case with nitrile/phenolic or vinyl butyral/phenolic adhesives, or impact-resistant molding compounds, additional silanes, particularly mercapto-functional silanes, have been found to impart greater bond strength than silanes that couple to the phenolic portion.



Polyurethane Coupling Reactions



Epoxy Coupling Reaction



Phenolic Coupling Reaction

Thermoplastics


Thermoplastics provide a greater challenge in promoting adhesion through silane coupling agents than thermosets. The silanes must react with the polymer and not the monomeric precursors, which not only limits avenues for coupling, but also presents additional problems in rheology and thermal properties during composite formulation. Moreover mechanical requirements here are stringently determined. Polymers that contain regular sites for covalent reactivity either in the backbone or in a pendant group include polydienes, polyvinylchloride, polyphenylene sulfide, acrylic homopolymers, maleic anhydride, acrylic, vinyl acetate, diene-containing copolymers, and halogen or chlorosulfonyl-modified homopolymers. A surprising number of these are coupled by aminoalkylsilanes. Chlorinated polymers readily form quaternary compounds while the carboxylate and sulfonate groups form amides and sulfonamides under process conditions. At elevated temperatures, the amines add across many double bonds although mercaptoalkylsilanes are the preferred coupling agents. The most widely used coupling agents, the aminoalkylsilanes, are not necessarily the best. Epoxysilanes, for example, are successfully used with acrylic acid and maleic acid copolymers.

Thermoplastic Condensation Polymers

The group of polymers that most closely approaches theoretical limits of composite strength does not appear to contain regular opportunities for covalent bond formation to substrate. Most of the condensation polymers including polyamides, polyesters, polycarbonates, and polysulfones are in this group. Adhesion is promoted by introducing high energy groups and hydrogen bond potential in the interphase area or by taking advantage of the relatively low molecular weight of these polymers, which results in a significant opportunity for end-group reactions. Aminoalkylsilanes, chloroalkylsilanes, and isocyanatosilanes are the usual candidates for coupling these resins. This group has the greatest mechanical strength of the thermoplastics, allowing them to replace the cast metals in such typical uses as gears, connectors and bobbins.




Scanning electron micrograph at a broken gear tooth from a non-coupled glass fiber/acetal composite. Note that cleavage occurred between fibers.



Scanning electron micrograph at a broken gear tooth from an aminosilane-coupled glass fiber/nylon 6/6 composite. Note how fibers have broken as well as matrix.



Chopped fiberglass strand sized with aminosilanes is a commonly used reinforcement for high temperature thermoplastics.




Thermoplastic Polyester Coupling Reaction

Polyolefins

The polyolefins and polyethers present no direct opportunity to covalent coupling. Until recently, the principal approach for composite formulation was to match the surface energy of the filler surface, by treating it with an alkyl-substituted silane, with that of the polymer. For optimum reinforcement, preferred resins should be of high molecular weight, linear, and have low melt viscosity. Approaches to improved composite strength have been through compatibility with long-chain alkylsilanes or aminosilanes. Far more effective is coupling with vinyl or methacryloxy groups, particularly if additional coupling sites are created in the resin by addition of peroxides. Dicumyl peroxide and bis(t-butylperoxy) compounds at levels of 0.15% to 0.25% have been introduced into polyethylene compounded with vinylsilane-treated glass fibers for structural composites or vinylsilane-treated clay for wire insulation. Increases of 50% in tensile and flexural properties have been observed in both cases when compared to the same silane systems without peroxides.

Another approach for coupling polypropylene and polyethylene is through silylsulfonylazides. Unlike azide bound to silicon, sulfonyl azides decompose above 150°C to form a molecule of nitrogen and a reactive nitrene that is capable of insertion into carbon-hydrogen bonds, forming sulfonamides, into carbon-carbon double bonds, forming triazoles, and into aromatic bonds, forming sulfonamides. Fillers are treated first with the silane and then the treated filler is fluxed rapidly with polymer melt.

Vinylsilanes are used in PE and EPDM insulated wire and cable



Polypropylene Coupling Reaction

Selecting a Silane Coupling Agent - Interphase Considerations

The space between homogeneous phases is sometimes called the interphase. In this region there is a steep gradient in local properties of the system. By treating a substrate with silanes the interphase can acquire specific surface energy, partition characteristics, mechanical and chemical properties.

Hydrophobicity and Wetting


Alkyl- and arylsilanes are not considered coupling agents. Surface modification with these non-functional materials can have profound effects on the interphase. They are used to alter surface energy or wetting characteristics of the substrate. In the simplest cases, methyltrichlorosilane, dimethyldichlorosilane, trimethylchlorosilane, their alkoxy derivatives, and hexamethyldisilazane are used to render substrates water repellent. For example, glassware can be dipped into a 5% to 10% solution of dimethyldiethoxysilane and heated for ten minutes at 120° C to render the surface hydrophobic. Laboratory pipettes and graduated cylinders so treated exhibit a flat meniscus and completely transfer aqueous solutions. GC packing of diatomaceous earth or silica are often treated with dimethyldichlorosilane or trimethylchlorosilane to reduce tailing. Masonry can be treated with propyl-, isobutyl- or octyltrialkoxysilanes to render it water repellent while glass surfaces of metal-glass capacitors treated with alkylsilanes exhibit reduced electrical leakage in humid conditions.

Silanes can alter the critical surface tension of a substrate in a well-defined manner. Critical surface tension is associated with the wettability or release qualities of a substrate. Liquids with a surface tension below the critical surface tension (γ_c) of a substrate will wet the surface, i.e., show a contact angle of 0 ($\cos\theta_c = 1$). The critical surface tension is unique for any solid, and is determined by plotting the cosine of the contact angles of liquids of different surface tensions and extrapolating to 1. The contact angle is given by Young's equation:

$$\gamma_{sv} - \gamma_{sl} = \cos\theta_e \gamma_{lv}$$

where γ_{sl} = interfacial surface tension, γ_{lv} = surface tension of liquid, and ($\gamma_{sv} = \gamma_l$ when $\gamma_{sl} = 0$ and $\cos \theta_e = 1$)

Contact Angle Defines Wettability



Critical surface tensions

	γ_c
Heptadecafluorodecyltrichlorosilane	12.0
Polytetrafluoroethylene	18.5
Methyltrimethoxysilane	22.5
Vinyltriethoxysilane	25
Paraffin wax	25.5
Ethyltrimethoxysilane	27.0
Propyltrimethoxysilane	28.5
Glass, soda-lime (wet)	30.0
Polychlorotrifluoroethylene	31.0
Polypropylene	31.0
Polyethylene	33.0
Trifluoropropyltrimethoxysilane	33.5
3-(2-aminoethyl)-aminopropyltrimethoxysilane	33.5
Polystyrene	34
Cyanoethyltrimethoxysilane	34
Aminopropyltriethoxysilane	35
Polyvinylchloride	39
Phenyltrimethoxysilane	40.0
Chloropropyltrimethoxysilane	40.5
Mercaptopropyltrimethoxysilane	41
Glycidoxypropyltrimethoxysilane	42.5
Polyethyleneterephthalate	43
Copper (dry)	44
Aluminum (dry)	45
Iron (dry)	46
Nylon 6/6	46
Glass, soda-lime (dry)	47
Silica, fused	78
Titanium dioxide (Anatase)	91
Ferric oxide	107
Tin oxide	111

Note: Critical surface tensions for silanes refer to treated surfaces.

Silane treatment has allowed control of thixotropic activity of silica and clays in grease and oil applications. In the reinforcement of thermosets and thermoplastics with glass fibers, one approach for optimizing reinforcement is to match the critical surface tension of the silylated glass surface to the surface tension of the polymer in its melt or uncured condition. This has been most helpful in resins with no obvious functionality such as polyethylene and polystyrene. Immobilization of cellular organelles, including mitochondria, chloroplasts, and microsomes, has been effected by treating silica with alkylsilanes of C_8 or greater substitution.

Chromatography

Octadecyl, cyanopropyl and branched tricocyl silanes provide bonded phases for liquid chromatography. Reverse-phase thin-layer chromatography can be accomplished by treating plates with dodecyltrichlorosilane.

By forming complexes of copper ions with amino-alkylsilylated surfaces, an interphase results that can selectively absorb ethylene, propylene and other gases.


Liquid Crystal Displays

The interphase can also impose orientation of the bulk phase. In liquid crystal displays, clarity and permanence of image are enhanced if the display can be oriented parallel or perpendicular to the substrate. The use of surfaces treated with octadecyl(3-(trimethoxysilyl)propyl) ammonium chloride (perpendicular) or methylaminopropyl-trimethoxysilane (parallel) has eliminated micromachining operations. The oriented crystalline domains often observed in reinforced nylons have also been attributed to orientation effects of the silane in the interphase.


Self-Assembled Monolayers (SAMs)

The perpendicular orientation of silanes with C_{10} or greater length can be utilized in micro-contact printing and other soft lithography methods. Here the silane may effect a simple differential adsorption, or if functionalized have a direct sensor effect.

Normal Phase HPLC of Carboxylic Acids with a C_{23} -Silane Bonded Phase




Orientation effects of silanes for passive LCDs



F. Kahn., Appl. Phys. Lett. 22, 386, 1973

Micro-Contact Printing Using SAMs



Special Topics

Dipodal Silanes

Functional dipodal silanes and combinations of non-functional dipodal silanes with functional silanes have significant impact on substrate bonding, hydrolytic stability and mechanical strength of many composites systems. They possess enabling activity in many coatings, particularly primer systems and aqueous immersion applications. The effect is thought to be a result of both the increased crosslink density of the interphase and a consequence of the fact that the resistance to hydrolysis of dipodal materials (with the ability to form six bonds to a substrate) is estimated at close to 100,000 times greater than conventional coupling agents (with the ability to form only three bonds to a substrate).


Both because dipodal silanes may not have functional groups identical to conventional coupling agents or because of economic considerations, conventional coupling agents are frequently used in combination with a non-functional dipodal silanes. In a typical application a dipodal material such as bis(triethoxysilyl)ethane (SIB1817.0) is combined at a 1:5 to 1:10 ratio with a traditional coupling agent. It is then processed in the same way as the traditional silane coupling agent.

Effect of dipodal -SiCH₂CH₂Si- on the bond strength of a crosslinkable ethylene-vinyl acetate primer formulation

Primer on metal 10% in <i>i</i> -PrOH	Wet adhesion to metals (N/cm)	
	Titanium	Cold-rolled steel
No silane	Nil	Nil
Methacryloxypropylsilane	0.25	7.0
Methacryloxypropylsilane + 10% dipodal	10.75	28.0 (cohesive failure)

90° peel strength after 2 h in 80°C water.


P. Pape et al, in *Silanes and Other Coupling Agents*, ed. K. Mittal, 1992, VSP, p105




Dipodal tetrasulfide silanes are used in “green” tires



Functional Dipodals




Non-Functional Dipodals




Linker Length

An important factor in controlling the effectiveness and properties of a coupled system is the linker between the organic functionality and the silicon atom. The linker length imposes a number of physical property and reactivity limitations. The desirability of maintaining the reactive centers close to the substrate are most important in sensor applications, in heterogeneous catalysis, fluorescent materials and composite systems in which the interfacing components are closely matched in modulus and coefficient of thermal expansion. On the other hand, inorganic surfaces can impose enormous steric constraints on the accessibility of organic functional groups in close proximity. If the linker length is long the functional group has greater mobility and can extend further from the inorganic substrate. This has important consequences if the functional group is expected to react with a single component in a multi-component organic or aqueous phases found in homogeneous and phase transfer catalysis, biological diagnostics or liquid chromatography. Extended linker length is also important in oriented applications such as self-assembled monolayers (SAMs). The typical linker length is three carbon atoms, a consequence of the fact that the propyl group is synthetically accessible and has good thermal stability.




Effect of linker length on the separation of aromatic hydrocarbons




T. Den et al, in "Silanes, Surfaces, Interfaces" D. Leyden ed., 1986 p403.


Silanes with short linker length



SIA0592.6



SIC2298.4



SII6453.8



SIP6723.7

Silanes with extended linker length




SIA0595.0




SIT8194.0



SIB1909.0



SIM6480.0




SIU9049.0



Cyclic Azasilanes

Volatile cyclic azasilanes are of particular interest in the surface modification of hydroxyl-containing surfaces, particularly inorganic surfaces such as nanoparticles and other nano-featured substrates. In these applications the formation of high functional density monolayers is critical. The cyclic azasilanes react with hydroxyl groups of a wide range of substrates at low temperatures by a ring-opening reaction that does not require water as a catalyst. Significantly, no byproducts of reaction form. The reactions of cyclic azasilanes are rapid at room temperature, even in the vapor phase. They also react rapidly at room temperature with isolated non-hydrogen bonded hydroxyls which do not undergo reaction with alkoxy-silanes under similar conditions. The three most common cyclic azasilanes structures are depicted. (see p.35)

Anhydrous deposition with Cyclic Azasilanes



Extent of reaction of organosilanes with fumed silica



M. Vedomuthu et al, J. Undergrad., Chem. Res., 1, 5, 2002


Thermal Stability of Silane Coupling Agents

The general order of thermal stability for silane coupling agents is depicted. Most commercial silane coupling agents have organic functionality separated from the silicon atom by three carbon atoms and are referred to as gamma-substituted silanes. The gamma-substituted silanes have sufficient thermal stability to withstand short-term process conditions of 350°C and long-term continuous exposure of 160°C. In some applications gamma-substituted silanes have insufficient thermal stability or other system requirements that can eliminate them from consideration. In this context, some comparative guidelines are provided for the thermal stability of silanes. Thermogravimetric Analysis (TGA) data for hydrolysates may be used for benchmarking. The specific substitution also plays a significant role in thermal stability. Electron withdrawing substitution reduces thermal stability, while electropositive groups enhance thermal stability.



Flexible multi-layer circuit boards for cell-phones utilize polyimide films coupled w/chloromethylaromatic silanes.

Relative Thermal Stability of Silanes



Thermal Stability of Silanes

SIA0025.0	$CH_3C(=O)CH_2CH_2Si(OC_2H_5)_3$	220°
SIC2271.0	$ClCH_2CH_2CH_2Si(OCH_3)_3$	360°
SIM6487.4	$H_2C=C(CH_3)C(=O)CH_2CH_2CH_2Si(OCH_3)_3$	395°
SIA0591.0	$H_2NCH_2CH_2NCH_2CH_2CH_2Si(OCH_3)_3$	390°
SIA0588.0	$H_2NCH_2CH_2NCH_2-\text{C}_6\text{H}_4-CH_2CH_2Si(OCH_3)_3$	435°
SIC2295.5	$ClCH_2-\text{C}_6\text{H}_4-CH_2CH_2Si(OCH_3)_3$	495°
SIA0599.1	$H_2N-\text{C}_6\text{H}_4-Si(OC_2H_5)_3$	485°
SIT8042.0	$CH_3-\text{C}_6\text{H}_4-Si(OCH_3)_3$	530°

25% weight loss of dried hydrolysates as determined by TGA


Aqueous Systems & Water-borne Silanes

Before most surface modification processes, alkoxy silanes are hydrolyzed forming silanol-containing species. The silanol-containing species are highly reactive intermediates which are responsible for bond formation with the substrate. In principal, if silanol species were stable, they would be preferred for surface treatments. Silanols condense with other silanols or with alkoxy silanes to form siloxanes. This can be observed when preparing aqueous treatment solutions. Initially, since most alkoxy silanes have poor solubility in water, two phases are observed. As the hydrolysis proceeds, a single clear phase containing reactive silanols forms. With aging, the silanols condense forming siloxanes and the solution becomes cloudy. Eventually, as molecular weight of the siloxanes increases, precipitation occurs.


Hydrolysis and condensation of alkoxy silanes is dependent on both pH and catalysts. The general objective in preparing aqueous solutions is to devise a system in which the rate of hydrolysis is substantially greater than the rate of condensation beyond the solubility limit of the siloxane oligomers. Other considerations are the work-time requirements for solutions and issues related to byproduct reactivity, toxicity or flammability. Stable aqueous solutions of silanes are more readily prepared if byproduct or additional alcohol is present in the solution since they contribute to an equilibrium condition favoring monomeric species.

Water-borne coupling agent solutions are usually free of VOCs and flammable alcohol byproducts. Most water-borne silanes can be described as hydroxyl-rich silsesquioxane copolymers. Apart from coupling, silane monomers are included to control water-solubility and extent of polymerization. Water-borne silanes act as primers for metals, additives for acrylic latex sealants and as coupling agents for siliceous surfaces.

Relative Hydrolysis Rates of Hydrolyzable Groups




Hydrolysis Profile of Phenylbis(2-methoxyethoxy)silanol




F. Osterholtz et al in Silanes and Other Coupling Agents ed K. Mittal, VSP, 1992, p119

Profile for Condensation of Silanols to Disiloxanes



E. Pohl et al in Silanes Surfaces and Interfaces ed., D. Leyden, Gordon and Breach, 1985, p481.



Water-borne Silsesquioxane Oligomers

Code	Functional Group	Mole %	Molecular Weight	Weight % in solution
WSA-7011	Aminopropyl	65-75	250-500	25-28
WSA-9911	Aminopropyl	100	270-550	22-25
WSA-7021	Aminoethylaminopropyl	65-75	370-650	25-28
WSAV-6511	Aminopropyl, Vinyl	60-65	250-500	25-28

Masked Silanes - Latent Functionality


Maximum bond strength in some adhesion and bonding systems requires that the organic functionality of a silane coupling agent becomes available during a discrete time period of substrate - matrix contact. Examples are epoxy adhesives in which reaction of the silane with the resin increases viscosity of an adhesive to the extent that substrate wet-out is inhibited and

pretreated fillers for composites which can react prematurely with moisture before melt compounding or vulcanization. A general approach is to mask the organic functionality of the silane which converts it to a storage-stable form and then to trigger the demasking with moisture, or heat concomitant with bonding or composite formation.


Masked Silanes - Moisture Triggered

Single-component liquid-cure epoxy adhesives and coatings employ dimethylbutylidene blocked amino silanes. These materials show excellent storage stability in resin systems, but are activated by moisture pro-

vided by water adsorbed on substrate surfaces or from humidity. Deblocking begins in minutes and is generally complete within two hours in sections with a diffusional thickness of less than 1mm.




Storage Stability of Epoxy Coating Solutions with blocked and unblocked aminosilanes




Epoxy Resin Solution: 50 parts bisphenol A epoxide, 5 parts SID4068.0 or SIA0610.0, 50 parts toluene.

Hydrolysis of Blocked Aminosilane

(SID4068.0/H₂O/THF = 1/10/20wt%)




An alternative is to use the moisture adsorbed onto fillers to liberate alcohol which, in turn, demasks the organic functionality.



Masked Silanes - Heat Triggered

Isocyanate functionality is frequently delivered to resin systems during elevated temperature bonding or melt processing steps. Demasking temperatures are typically 160-200°C.



PLEASE INQUIRE ABOUT BULK QUANTITIES

Coupling Agents for Metal Substrates

The optimum performance of silane coupling agents is associated with siliceous substrates. While the use of silanes has been extended to metal substrates, both the effectiveness and strategies for bonding to these less-reactive substrates vary. Four approaches of bonding to metals have been used with differing degrees of success. In all cases, selecting a dipodal or polymeric silane is preferable to a conventional trialkoxy silane.


Metals that form hydrolytically stable surface oxides, e.g. aluminum, tin, titanium. These oxidized surfaces tend to have sufficient hydroxyl functionality to allow coupling under the same conditions applied to the siliceous substrates discussed earlier.

Metals that form hydrolytically or mechanically unstable surface oxides, e.g. iron, copper, zinc. These oxidized surfaces tend to dissolve in water leading to progressive corrosion of the substrate or form a passivating oxide layer without mechanical strength. The successful strategies for coupling to these substrates typically involves two or more silanes. One silane is a chelating agent such as a diamine, polyamine or polycarboxylic acid. A second silane is selected which has a reactivity with the organic component and reacts with the first silane by co-condensation. If a functional dipodal or polymeric silane is not selected, 10-20% of a non-functional dipodal silane typically improves bond strength.

Metals that do not readily form oxides, e.g. nickel, gold and other precious metals. Bonding to these substrates requires coordinative bonding, typically a phosphine, sulfur (mercapto), or amine functional silane. A second silane is selected which has a reactivity with the organic component. If a functional dipodal or polymeric silane is not selected, 10-20% of a non-functional dipodal silane typically improves bond strength.




SIP6926.2



Octylsilane adsorbed on titanium

figure courtesy of M. Banaszak-Holl

Metals that form stable hydrides, e.g. titanium, zirconium, nickel. In a significant departure from traditional silane coupling agent chemistry, the ability of certain metals to form so-called amorphous alloys with hydrogen is exploited in an analogous chemistry in which hydride functional silanes adsorb and then coordinate with the surface of the metal. Most silanes of this class possess only simple hydrocarbon substitution such as octylsilane. However they do offer organic compatibility and serve to markedly change wet-out of the substrate. Both hydride functional silanes and treated metal substrates will liberate hydrogen in the presence of base or with certain precious metals such as platinum and associated precautions must be taken. (see p 53.)



Coupling Agents for Metals*

Metal	Class	Screening Candidates	
Copper	Amine	SSP-060	SIT8398.0
Gold	Sulfur	SIT7908.0	SIP6926.2
	Phosphorus	SID4558.0	SIB1091.0
Iron	Amine	SIB1834.0	WSA-7011
	Sulfur	SIB1824.6	SIM6476.0
Tin	Amine	SIB1835.5	
Titanium	Epoxy	SIG5840.0	SIE6668.0
	Hydride	SIU9048.0	
Zinc	Amine	SSP-060	SIT8398.0
	Carboxylate	SIT8402.0	SIT8192.6

*These coupling agents are almost always used in conjunction with a second silane with organic reactivity or a dipodal silane.


Difficult Substrates

Silane coupling agents are generally recommended for applications in which an inorganic surface has hydroxyl groups and the hydroxyl groups can be converted to stable oxane bonds by reaction with the silane. Substrates such as calcium carbonate, copper and ferrous alloys, and high phosphate and sodium glasses are not recommended substrates for silane coupling agents. In cases where a more appropriate technology is not available a number of strategies have been devised which exploit the organic functionality, film-forming and crosslinking properties of silane coupling agents as the primary mechanism for substrate bonding in place of bonding through the silicon atom. These approaches frequently involve two or more coupling agents.

Calcium carbonate fillers and marble substrates do not form stable bonds with silane coupling agents. Applications of mixed silane systems containing a dipodal silane or tetraethoxysilane in combination with an organofunctional silane frequently increases adhesion. The adhesive mechanism is thought to be due to the low molecular weight and low surface energy of the silanes which allows them initially to spread to thin films and penetrate porous structures followed by the crosslinking which results in the formation of a silica-rich encapsulating network. The silica-rich encapsulating network is then susceptible to coupling chemistry comparable to siliceous substrates. Marble and calciferous substrates can also benefit from the inclusion of anhydride-functional silanes which, under reaction conditions, form dicarboxylates that can form salts with calcium ions.

Metals and many metal oxides can strongly adsorb silanes if a chelating functionality such as diamine or dicarboxylate is present. A second organofunctional silane with reactivity appropriate to the organic component must be present. Precious metals such as gold and rhodium form weak coordination bonds with phosphine and mercaptan functional silanes.

High phosphate and sodium content glasses are frequently the most frustrating substrates. The primary inorganic constituent is silica and would be expected to react readily with silane coupling agents. However alkali metals and phosphates not only do not form hydrolytically stable bonds with silicon, but, even worse, catalyze the rupture and redistribution of silicon-oxygen bonds. The first step in coupling with these substrates is the removal of ions from the surface by extraction with deionized water. Hydrophobic dipodal or multipodal silanes are usually used in combination with organofunctional silanes. In some cases polymeric silanes with multiple sites for interaction with the substrate are used. Some of these, such as the polyethylenimine functional silanes can couple to high sodium glasses in an aqueous environment.



Substrates with low concentrations of non-hydrogen bonded hydroxyl groups, high concentrations of calcium, alkali metals or phosphates pose challenges for silane coupling agents.

Removing Surface Impurities

Eliminating non-bonding metal ions such as sodium, potassium and calcium from the surface of substrates can be critical for stable bonds. Substrate selection can be essential. Colloidal silicas derived from tetraethoxysilane or ammonia sols perform far better than those derived from sodium sols. Bulk glass tends to concentrate impurities on the surface during fabrication. Although sodium concentrations derived from bulk analysis may seem acceptable, the surface concentration is frequently orders of magnitude higher. Surface impurities may be reduced by immersion in 5% hydrochloric acid for 4 hours, followed by a deionized water rinse, and then immersion in deionized water overnight followed by drying.

Oxides with high isoelectric points can adsorb carbon dioxide, forming carbonates. These can usually be removed by a high temperature vacuum bake.

Increasing Hydroxyl Concentration

Hydroxyl functionalization of bulk silica and glass may be increased by immersion in a 1:1 mixture of 50% aqueous sulfuric acid : 30% hydrogen peroxide for 30 minutes followed by rinses in D.I. water and methanol and then air drying. Alternately, if sodium ion contamination is not critical, boiling with 5% aqueous sodium peroxodisulfate followed by acetone rinse is recommended¹.
1. K. Shirai et al, J. Biomed. Mater. Res. 53, 204, 2000.

Catalyzing Reactions in Water-Free Environments

Hydroxyl groups without hydrogen bonding react slowly with methoxy silanes at room temperature. Ethoxy silanes are essentially non-reactive. The methods for enhancing reactivity include transesterification catalysts and agents which increase the acidity of hydroxyl groups on the substrate by hydrogen bonding. Transesterification catalysts include tin compounds such as dibutyltin diacetoxystyrene and titanates such as titanium isopropoxide. Incorporation of transesterification catalysts at 2-3 weight % of the silane effectively promotes reaction and deposition in many instances. Alternatively, amines can be premixed with solvents at 0.01-0.5 weight % based on substrate prior or concurrent to silane addition. Volatile primary amines such as butylamine can be used, but are not as effective as tertiary amines such as benzyltrimethylamine or diamines such as ethylenediamine. The more effective amines, however, are more difficult to remove after reaction¹.
1. S. Kanan et al, Langmuir, 18, 6623, 2002.

Hydroxylation by Water Plasma & Steam Oxidation

Various metals and metal oxides including silicon and silicon dioxide can achieve high surface concentrations of hydroxyl groups after exposure to H₂O/O₂ in high energy environments including steam at 1050° and water plasma¹.

1. N. Alcanter et al, in "Fundamental & Applied Aspects of Chemically Modified Surfaces" ed. J. Blitz et al, 1999, Roy. Soc. Chem., p212.

Applying Silanes

Deposition from aqueous alcohol solutions is the most facile method for preparing silylated surfaces. A 95% ethanol-5% water solution is adjusted to pH 4.5-5.5 with acetic acid. Silane is added with stirring to yield a 2% final concentration. Five minutes should be allowed for hydrolysis and silanol formation. Large objects, e.g. glass plates, are dipped into the solution, agitated gently, and removed after 1-2 minutes. They are rinsed free of excess materials by dipping briefly in ethanol. Particles, e.g. fillers and supports, are silylated by stirring them in solution for 2-3 minutes and then decanting the solution. The particles are usually rinsed twice briefly with ethanol. Cure of the silane layer is for 5-10 mins at 110°C or 24 hours at room temperature (<60% relative humidity).

Fig. 1 Reactor for slurry treatment of powders. Separate filtration and drying steps are required.



Deposition from aqueous solution is employed for most commercial fiberglass systems. The alkoxysilane is dissolved at 0.5-2.0% concentration in water. For less soluble silanes, 0.1% of a non-ionic surfactant is added prior to the silane and an emulsion rather than a solution is prepared. The solution is adjusted to pH 5.5 with acetic acid. The solution is either sprayed onto the substrate or employed as a dip bath. Cure is at 110-120°C for 20-30 minutes. Stability of aqueous silane solutions varies from 2-12 hours for the simple alkyl silanes. Poor solubility parameters limit the use of long chain alkyl and aromatic silanes by this method. Distilled water is not necessary, but water containing fluoride ions must be avoided.

Bulk deposition onto powders, e.g. filler treatment, is usually accomplished by a spray-on method. It assumes that the total amount of silane necessary is known and that sufficient adsorbed moisture is present on the filler to cause hydrolysis of the silane. The silane is prepared as a 25% solution in alcohol. The powder is placed in a high intensity solid mixer, e.g. twin cone mixer with intensifier. The methods are most effective. If the filler is dried in trays, care must be taken to avoid wicking or skinning of the top layer of treated material by adjusting heat and air flow.

Fig. 2 Vacuum tumble dryers can be used for slurry treatment of powders.



Integral blend methods are used in composite formulations. In this method the silane is used as a simple additive. Composites can be prepared by the addition of alkoxysilanes to dry-blends of polymer and filler prior to compounding. Generally 0.2 to 1.0 weight percent of silane (of the total mix) is dispersed by spraying the silane in an alcohol carrier onto a preblend. The addition of the silane to non-dispersed filler is not desirable in this technique since it can lead to agglomeration. The mix is dry-blended briefly and then melt compounded. Vacuum devolatilization of byproducts of silane reaction during melt compounding is necessary to achieve optimum properties. Properties are sometimes enhanced by adding 0.5-1.0% of tetrabutyl titanate or benzyldimethylamine to the silane prior to dispersal.

Anhydrous liquid phase deposition of chlorosilanes, methoxysilanes, aminosilanes and cyclic azasilanes is preferred for small particles and nano-featured substrates. Toluene, tetrahydrofuran or hydrocarbon solutions are prepared containing 5% silane. The mixture is refluxed for 12-24 hours with the substrate to be treated. It is washed with the solvent. The solvent is then removed by air or explosion-proof oven drying. No further cure is necessary. This reaction involves a direct nucleophilic displacement of the silane chlorines by the surface silanol. If monolayer deposition is desired, substrates should be predried at 150°C for 4 hours. Bulk deposition results if adsorbed water is present on the substrate. This method is cumbersome for large scale preparations and rigorous controls must be established to ensure reproducible results. More reproducible coverage is obtained with monochlorosilanes.

Chlorosilanes can also be deposited from alcohol solution. Anhydrous alcohols, particularly ethanol or isopropanol are preferred. The chlorosilane is added to the alcohol to yield a 2-5% solution. The chlorosilane reacts with the alcohol producing an alkoxysilane and HCl. Progress of the reaction is observed by halt of HCl evolution. Mild warming of the solution (30-40°C) promotes completion of the reaction. Part of the HCl reacts with the alcohol to produce small quantities of alkyl halide and water. The water causes formation of silanols from alkoxysilanes. The silanols condense on the substrate. Treated substrates are cured for 5-10 mins. at 110°C or allowed to stand 24 hours at room temperature.



Fig. 3 Twin-cone blenders with intensive mixing bars are used for bulk deposition of silanes onto powders.

Applying Silanes

Vapor Phase Deposition

Silanes can be applied to substrates under dry aprotic conditions by chemical vapor deposition methods. These methods favor monolayer deposition. Although under proper conditions almost all silanes can be applied to substrates in the vapor phase, those with vapor pressures >5 torr at 100°C have achieved the greatest number of commercial applications. In closed chamber designs, substrates are supported above or adjacent to a silane reservoir and the reservoir is heated to sufficient temperature to achieve 5mm vapor pressure. Alternatively, vacuum can be applied until silane evaporation is observed. In still another variation the silane can be prepared as a solution in toluene, and the toluene brought to reflux allowing sufficient silane to enter the vapor phase through partial pressure contribution. In general, substrate temperature should be maintained above 50° and below 120° to promote reaction. Cyclic azasilanes deposit the quickest—usually less than 5 minutes. Amine functional silanes usually deposit rapidly (within 30 minutes) without a catalyst. The reaction of other silanes requires extended reaction times, usually 4-24 hours. The reaction can be promoted by addition of catalytic amounts of amines.

Spin-On

Spin-On applications can be made under hydrolytic conditions which favor maximum functionalization and polylayer deposition or dry conditions which favor monolayer deposition. For hydrolytic deposition 2-5% solutions are prepared (see deposition from aqueous alcohol). Spin speed is low, typically 500 rpm. Following spin-deposition a hold period of 3-15 minutes is required before rinse solvent. Dry deposition employs solvent solutions such as methoxypropanol or ethyleneglycol monoacetate (EGMA). Aprotic systems utilize toluene or THF. Silane solutions are applied at low speed under a nitrogen purge. If strict monolayer deposition is preferred, the substrate should be heated to 50° . In some protocols, limited polylayer formation is induced by spinning under an atmospheric ambient with 55% relative humidity.

Spray application

Formulations for spray applications vary widely depending on end-use. They involve alcohol solutions and continuously hydrolyzed aqueous solutions employed in architectural and masonry applications. The continuous hydrolysis is effected by feeding mixtures of silane containing an acid catalyst such as acetic acid into a water stream by means of a venturi (aspirator). Stable aqueous solutions (see water-borne silanes), mixtures of silanes with limited stability (4-8 hours) and emulsions are utilized in textile and fiberglass applications. Complex mixtures with polyvinyl acetates or polyesters enter into the latter applications as sizing formulations.

Figure 4.
Apparatus for vapor phase silylation.





Figure 5.
Spin-coater for deposition on wafers.



Figure 6.
Spray application of silanes on large structures.
















Figure 7.
Spray & contact roller application of silanes on fiberglass.





Acrylate-silanes in dental restorative composites.

Silane Coupling Agents for Thermosets Selection Chart

	Coupling Agent Class	Suggestions for Primary Screening
Acrylate, UV cure 	Acrylate Vinyl/Olefin	SIA0200.0 SIM6487.4 SIS6964.0
Diallylphthalate 	Amine Vinyl/Olefin	SIA0591.0 SIA0610.0 SIS6964.0
Epoxy 	Amine Anhydride Epoxy	SIA0591.0 SIT8398.0 SIT8192.6 SIG5840.0
Epoxy, UV Cure 	Amine Epoxy	SIA0591.0 SIT8398.0 SIE4668.0 SIE4670.0
Polyimide 	Amine Halogen Dipodal	SIA0599.2 SIA0591.0 SIC2295.5 SIC2296.2 SIB1833.0
Furan 	Amine Epoxy	SIA0611.0 SIA0599.0 SIG5840.0
Melamine 	Amine Hydroxyl Dipodal	SIA0611.0 SIA0599.0 SIB1140.0 SIB1833.0 SIT8717.0
Parylene 	Halogen Vinyl/Olefinic Dipodal	SIC2295.5 SIS6990.0 SIM6487.4 SIB1832.0 VMM-010
Phenol-formaldehyde 	Amine Epoxy	SIA0611.0 SIT8187.5 SIE4670.0 SIG5840.0
Methylmethacrylate, cast 	Acrylate Amine	SIM6487.4 SIA0200.0 SIB1828.0
Polyester, unsaturated 	Acrylate Vinyl/Olefin	SIM6487.4 SIS6994.0 SIV9112.0
Urea-formaldehyde 	Amine Hydroxyl	SIA0610.0 SIU9055.0 SIB1140.0
Urethane 	Amine Isocyanate Sulfur	SIA0610.0 SIM6500.0 SII6455.0 SIM6476.0

Diamine-silanes
couple
polycarbonate
in CDs



Silane Coupling Agents for Thermoplastics Selection Chart


	Coupling Agent Class	Suggestions for Primary Screening
Polyacetal	$\left[-\text{CH}_2\text{O}- \right]_n$ Vinyl/Olefin	SIS6994.0
Polyacrylate	$\left[\begin{array}{c} \text{O} \\ \parallel \\ \text{C}-\text{OCH}_3 \\ \\ \text{CH}_2-\text{C} \\ \\ \text{H} \end{array} \right]_n$ Amine	SIU9058.0 SIA0610.0
Polyamide	$\left[-\text{NH}(\text{CH}_2)_m\text{C}(=\text{O})- \right]_n$ Amine Dipodal Water-borne	SIA0610.0 SIB1834.1 WSA-7011 SIA0614.0 SSP-060
Polyamide-imide	$\left[\begin{array}{c} \text{O} \\ \parallel \\ \text{N} \\ \\ \text{C} \\ \\ \text{H} \end{array} \right]_n$ Amine Halogen	SIA0610.0 SIC2295.5
Polybutylene terephthalate	$\left[\begin{array}{c} \text{O} \\ \parallel \\ \text{C} \\ \\ \text{C}_6\text{H}_4 \\ \\ \text{C} \\ \parallel \\ \text{O} \end{array} \right]_n$ Amine Isocyanate	SIA0610.0 SII6455.0
Polycarbonate	$\left[\begin{array}{c} \text{O} \\ \parallel \\ \text{C} \\ \\ \text{C}_6\text{H}_4 \\ \\ \text{C} \\ \parallel \\ \text{O} \end{array} \right]_n$ Amine	SIA0591.0 SIA0610.0
Polyether ketone	$\left[\begin{array}{c} \text{O} \\ \parallel \\ \text{C} \\ \\ \text{C}_6\text{H}_4 \\ \\ \text{C}_6\text{H}_4 \\ \\ \text{O} \end{array} \right]_n$ Amine Dipodal	SIA0591.0 SIT8717.0
Polyethylene	$\left[-\text{CH}_2\text{CH}_2- \right]_n$ Amine Vinyl/Olefin	SIA0591.0 SSP-055 SIT8398.0 SIV9112.0
Polyphenylene sulfide	$\left[\begin{array}{c} \text{C}_6\text{H}_4 \\ \\ \text{S} \\ \\ \text{C}_6\text{H}_4 \end{array} \right]_n$ Amine Halogen Sulfur	SIA0605.0 SIC2295.5 SIM6476.0
Polypropylene	$\left[\begin{array}{c} \text{CH}_3 \\ \\ -\text{CH}_2\text{CH}- \end{array} \right]_n$ Acrylate Azide Vinyl/Olefin	SIM6487.4 SIA0780.0 VEE-005 SSP-055
Polystyrene	$\left[\begin{array}{c} \text{CH}_2\text{CH} \\ \\ \text{C}_6\text{H}_5 \end{array} \right]_n$ Acrylate Dipodal	SIM6487.4 SIB1831.0
Polysulfone	$\left[\begin{array}{c} \text{O} \\ \parallel \\ \text{C} \\ \\ \text{C}_6\text{H}_4 \\ \\ \text{C} \\ \parallel \\ \text{O} \end{array} \right]_n$ Amine	SIA0591.0 SIU9055.0
Polyvinyl butyral	$\left[\begin{array}{c} \text{CH}_2 \\ \\ \text{C} \\ \\ \text{O} \\ \\ \text{CH}_2\text{CH}_2\text{CH}_3 \end{array} \right]_n$ Amine	SIA0611.0 SIU9058.0
Polyvinyl chloride	$\left[\begin{array}{c} \text{Cl} \\ \\ -\text{CH}_2\text{CH}- \end{array} \right]_n$ Amine Sulfur	SIA0605.0 SIM6474.0 SIB1825.0

Silane Coupling Agents for Sealants & Elastomers Selection Chart

Water-borne aminosilanes increase bonding of acrylic latex sealants



Acrylic latex



Coupling Agent Class






Suggestions for Primary Screening

		Acrylate Vinyl/Olefin Water-borne	SIM6487.4 SIV9210.0 WSA-7021	SIV9218.0 WSA-6511
Butyl	$\left[\text{CH}_2\text{CH}=\text{CHCH}_2 \right]_n$	Acrylate Sulfur Vinyl/Olefin	SIM6487.4 SIB1825.0 SSP-055	SIM6476.0 VEE-005
Epichlorohydrin	$\left[\begin{array}{c} \text{OCH}_2\text{CH} \\ \\ \text{CH}_2\text{Cl} \end{array} \right]_n$	Amine Sulfur	SIA0605.0 SIM6474.0	
Fluorocarbon	$\text{---}(\text{CF}_2\text{CF}_2)_m(\text{CH}_2\text{CF}_2)_n\text{---}$	Amine Dipodal	SIB1834.1 SIT8717.0	
Isoprene	$\left[\begin{array}{c} \text{CH}_3 \\ \\ \text{CH}_2\text{C}=\text{CHCH}_2 \end{array} \right]_n$	Sulfur Vinyl/Olefin	SIM6474.0 SSP-055	SIM6476.0 VEE-005
Neoprene	$\left[\begin{array}{c} \text{Cl} \\ \\ \text{CH}_2\text{C}=\text{CHCH}_2 \end{array} \right]_n$	Sulfur Vinyl/Olefin	SIM6474.0 SSP-055	SIM6476.0 VEE-005
Nitrile	$\left[\begin{array}{c} \text{CN} \\ \\ \text{CH}_2\text{CH}-\text{CH}_2-\text{CH}=\text{CH} \end{array} \right]_n$	Epoxy Sulfur	SIG5840.0 SIB1825.0	
Polysulfide	$\left[\text{CH}_2\text{CH}_2\text{S} \right]_n$	Epoxy Sulfur	SIG5840.0 SIB1825.0	SIM6476.0
SBR	$\left[\begin{array}{c} \text{---} \\ \\ \text{CH}_2\text{CH}-\text{CH}_2-\text{CH}=\text{CH} \\ \\ \text{C}_6\text{H}_5 \end{array} \right]_n$	Amine Sulfur	SIA0605.0 SIB1825.0	SIM6486.0
Silicone (hydroxyl terminated)	$\text{HO}-\text{Si}(\text{CH}_3)_2-\text{O}-\left(\text{Si}(\text{CH}_3)_2-\text{O} \right)_n-\text{Si}(\text{CH}_3)_2-\text{OH}$	Amine Vinyl/Olefin Dipodal	SIA0605.0 SIV9098.0 SIB1824.0	SIA0589.0 VMM-010
Silicone (vinyl terminated)	$\text{H}_2\text{C}=\text{CH}-\text{Si}(\text{CH}_3)_2-\text{O}-\left(\text{Si}(\text{CH}_3)_2-\text{O} \right)_n-\text{Si}(\text{CH}_3)_2-\text{CH}=\text{CH}_2$	Acrylate Vinyl/Olefin	SIM6487.4 SIA0540.0	VMM-010



aldehyde-, amino-, and hydroxyl-silanes couple DNA in array technology

Silane Coupling Agents for Biomaterials Selection Chart

	<i>Site/Type</i>	<i>Coupling Class</i>	<i>Co-reactant</i>	<i>Suggestions for Screening</i>	
	<i>Oligonucleotides</i>	<i>hydroxyl diamine</i>	<i>cobalt ethylenediamine</i>	<i>SIB1140.0 SIA0591.0</i>	
				<i>G. McGall et al, J. Am. Chem. Soc., 119, 5081, 1997. F. Chow, in "Silylated Surfaces" D. Leyden ed., Gordon & Breach, 1978, p.301.</i>	
	<i>DNA</i>	<i>terminal favored pendant amine pendant amine pendant amine</i>	<i>vinyl/olefin aldehyde diamine epoxy</i>	<i>SIO6708.0 SIT8194.0 SIA0594.0 SIE4675.0</i> <i>SIU9049.0 SID3543.0 SIG5838.0</i>	
				<i>A. Bensimon, Science, 265, 2096, 1994. J. Grobe et al, J. Chem. Soc. Chem. Commun, 2323, 1995. C. Kneuer et al, Int'l J. Pharmaceutics, 196(2), 257, 2000.</i>	
		<i>Protein</i>	<i>lysine lysine lysine cysteine tyrosine heparinated immunoglobulin antibody</i>	<i>aldehyde amine amine sulfur nitrobenzamide amine/quat pyridyl-thio cyano</i>	<i>SIT8194.0 SIA0611.0 SIA0611.0 SIM6476.0 SIT8191.0 SSP-060 SIP6926.4 SIC2456.0</i> <i>SIA0595.0 SIA0599.0 SIT8415.0</i>
					<i>J. Grobe et al, J. Chem. Soc. Chem. Commun, 2323, 1995. H. Weetall, US Pat. 3,652,761. G. Royer, CHEMTECH, 4, 699, 1974. S. Bhatia et al, Anal. Biochem., 178, 408, 1989. J. Venter et al, Proc. Nat. Acad. Soc., 69(5), 1141, 1972. R. Merker et al, Proc. Artificial Heart Prog. Conf., June 9-13, 1969 HEWNIH, p29. S. Falipou, Fundamental & Applied Aspects of Chemically Modified Surfaces, p389, 1999.</i>
		<i>Cell-Organelle</i>	<i>chloroplast mitochondria</i>	<i>alkyl alkyl</i>	<i>SIO6645.0 SIO6645.0</i>
		<i>mitochondria on silica bead</i>			<i>B. Arkles et al, in "Silylated Surfaces" D. Leyden ed., Gordon & Breach, 1978, p363. B. Arkles et al, J. Biol. Chem., 250, 8856, 1975.</i>
		<i>Whole Cell</i>	<i>erythrocytes</i>	<i>short alkyl</i>	<i>SIE4901.4</i>
		<i>erythrocytes on glass wall</i>			<i>B. Arkles et al, in "Silylated Surfaces" D. Leyden ed., Gordon & Breach, 1978, p363.</i>
	<i>Whole Cell (causing lysis)</i>	<i>procaryotic</i>	<i>alkyl-quat</i>	<i>SIO6620.0 SID3392.0</i>	
				<i>W. White et al in "Silanes, Surfaces & Interfaces" ed. D. Leyden, Gordon & Breach, 1986, p. 107.</i>	
	<i>Tissue</i>	<i>histological samples</i>		<i>SIA0611.0 SIA0610.0</i>	



SILANE COUPLING AGENT PROPERTIES



Lactase is immobilized with aminosilanes and glutaraldehyde.



Epoxy-silanes are essential for performance of epoxy resin encapsulants for microchips.



Methacrylate-silanes couple fiberglass to unsaturated polyester in corrosion resistant rooftop ductwork at Gelest, Inc.

Acrylate & Methacrylate functional	26
Aldehyde functional	27
Amino functional	28
Anhydride functional.	36
Azide functional	36
Carboxylate, Phosphonate and Sulfonate functional	36
Epoxy functional.	37
Ester functional.	38
Halogen functional	38
Hydroxyl functional.	40
Isocyanate and Masked Isocyanate functional.	41
Phosphine and Phosphate functional	42
Sulfur functional.	43
Vinyl and Olefin functional	45
Multi-functional and Polymeric Silanes.	49
Water-borne Coupling Agents	49
Non-functional Dipodal Silanes.	50
UV Active and Fluorescent Silanes.	51
Chiral Silanes	52
Biomolecular Probes	53
Silyl Hydrides	53

Commercial Status - produced on a regular basis for inventory

Developmental Status - available to support development and commercialization




Methacrylate-silanes
couple fiberglass
to unsaturated
polyester


Acrylate & Methacrylate Functional Silanes

name MW bp/mm (mp) D₄²⁰ n_D²⁰


Acrylate & Methacrylate Functional Silanes - Trialkoxy




SIA0200.0
(3-ACRYLOXYPROPYL)TRIMETHOXY-
SILANE, 95% inhibited with MEHQ
C₉H₁₈O₅Si
aqueous solutions more stable than methacrylate analog
coupling agent for epoxies, UV cure coatings; employed in optical fiber coatings¹.
1. M. Yokoshima et al, CA113, 15746d; Jap. Pat. 02133338, 1990
[4369-14-6] TSCA-S HMIS: 3-1-1-X store <5° 25g/\$48.00 100g/\$156.00 2.0kg/\$580.00




SIM6487.4
METHACRYLOXYPROPYLTRIMETHOXY-
SILANE MEMO inhibited with MEHQ, HQ
C₁₀H₂₀O₅Si
viscosity: 2 cSt.
copolymerization parameters-e,Q: 0.07, 2.7
widely used coupling agent for unsaturated polyester-fiberglass composites¹.
copolymerized with styrene in formation of sol-gel composites².
1. B. Arkles, Chemtech, 7, 713, 1977
2. Y. Wei et al, J. Mater. Res., 8, 1143, 1993
[2530-85-0] TSCA HMIS: 3-2-1-X store <5° 100g/\$10.00 2.0kg/\$124.00 18kg/\$630.00




SIA0180.0
N-(3-ACRYLOXY-2-HYDROXYPROPYL)-3-
AMINOPROPYLTRIETHOXYSILANE, 50% in ethanol
C₁₅H₃₁NO₆Si inhibited with MEHQ
flashpoint: 8°C (48°F)
[123198-57-2] HMIS: 3-4-1-X store <5° 25g/\$140.00




SIM6480.8
O-(METHACRYLOXYETHYL)-N-(TRIETHOXY-
SILYLPROPYL)URETHANE, 90%
C₁₆H₃₁NO₇Si inhibited with MEHQ
flashpoint: 8°C (48°F)
[115396-93-5] HMIS: 3-2-1-X store <5° 25g/\$42.00 100g/\$136.00




SIM6481.1
N-(3-METHACRYLOXY-2-HYDROXYPROPYL)-3-
AMINOPROPYLTRIETHOXYSILANE, 50% in ethanol
C₁₆H₃₃NO₆Si inhibited with MEHQ
flashpoint: 8°C (48°F)
employed in conservation/consolidation of stone¹.
1. G. Wheeler, in "Ninth Int'l Cong. on Deteriorat'n and Conservat'n of Stone" ed. V Fassina, 2, 541, Elsevier 2000.
[96132-98-8] HMIS: 3-4-1-X store <5° 25g/\$49.00 100g/\$160.00



SIM6482.0
METHACRYLOXYMETHYLTRIETHOXYSILANE
C₁₁H₂₂O₅Si inhibited with MEHQ
treatment of fumed silica in acrylic casting compositions accelerates polymerization¹.
1. E. Morozova et al, CA 95,98753g; Plast. Massy, 7, 1981
[5577-72-0] HMIS: 3-2-1-X store <5° 10g/\$40.00 50g/\$160.00













SIM6483.0
METHACRYLOXYMETHYLTRIMETHOXY-
SILANE
C₈H₁₆O₅Si inhibited with MEHQ
modification of novolac resins afford bilevel resists having attributes of trilevel resists¹.
1. E. Reichmanis et al, US Pat. 4,481,049,1984
[54586-78-6] HMIS: 3-2-1-X store <5° 10g/\$32.00 50g/\$128.00



SIM6487.3
METHACRYLOXYPROPYLTRIETHOXYSILANE
C₁₃H₂₆O₅Si inhibited with MEHQ
flashpoint: 128°C (262°F)
[21142-29-0] HMIS: 3-1-1-X store <5° 10g/\$39.00 50g/\$156.00

Commercial

Developmental

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
Acrylate & Methacrylate Functional Silanes - Dialkoxy					
	SIA0198.0 (3-ACRYLOXYPROPYL)METHYLDIMETHOXY SILANE, 95% inhibited w/ MEHQ C ₉ H ₁₈ O ₄ Si employed in fabrication of photoimageable, low shrinkage multimode waveguides ¹ . 1. C. Xu et al, Chem. Mater., 8, 2701, 1996	218.33	65°/0.35	1.0	1.431
	[13732-00-8] HMIS: 3-2-1-X store <5°	50g/\$46.00	250g/\$184.00		
	SIM6481.43 (METHACRYLOXYMETHYL)METHYL- DIETHOXY SILANE C ₁₀ H ₂₀ O ₄ Si [121177-93-3] HMIS: 2-2-1-X store <5°	232.4	221°	0.977	
	flashpoint: 88°C (190°F) TOXICITY oral-rat, LD50: >2000mg/kg	25g/\$26.00	100g/\$84.00		
	SIM6481.46 (METHACRYLOXYMETHYL)METHYL- DIMETHOXY SILANE C ₈ H ₁₆ O ₄ Si viscosity: 1.4 cSt [3978-58-3] HMIS: 3-2-1-X store <5°	204.30	205°	1.020	1.4274
	flashpoint: 82°C (180°F) autoignition temp: 300°C TOXICITY oral-rat, LD50: >2000mg/kg	25g/\$22.00	100g/\$72.00		
	SIM6486.8 METHACRYLOXYPROPYLMETHYLDI- ETHOXY SILANE, 95% inhibited w/ MEHQ C ₁₂ H ₂₄ O ₄ Si [65100-04-1] HMIS: 3-1-1-X store <5°	260.40	95°/1	0.965	1.433
	flashpoint 136°C (277°F)	10g/\$39.00	50g/\$156.00		
	SIM6486.9 METHACRYLOXYPROPYLMETHYLDI- METHOXY SILANE, 95% inhibited w/MEHQ C ₁₀ H ₂₀ O ₄ Si monomer for hybrid inorganic-organic composites ¹ . 1. R. Taylor-Smith, Polym. Mat. Sci. Eng., Preprints, 77, 503, 1997	235.69	83°/3	1.00	1.4351
	flashpoint: 115°C (190°F)	[14513-34-9] HMIS: 3-2-1-X store <5°	25g/\$45.00	100g/\$146.00	
Acrylate & Methacrylate Functional Silanes - Monoalkoxy					
	SIM6486.4 METHACRYLOXYPROPYLDIMETHYLETHOXY- SILANE, 95% inhibited with MEHQ C ₁₁ H ₂₂ O ₃ Si [13731-98-1] HMIS: 3-2-1-X store <5°	230.38	75-6°/0.4	0.926	1.4371
	10g/\$78.00				
	SIM6486.5 METHACRYLOXYPROPYLDIMETHYL- METHOXY SILANE, 95% inhibited with MEHQ C ₁₀ H ₂₀ O ₃ Si [66753-64-8] HMIS: 3-2-1-X store <5°	216.35	70-2°/0.5	0.944	1.4381
	10g/\$48.00	50g/\$192.00			
Aldehyde Functional Silanes					
Aldehyde Functional Silanes - Trialkoxy					
	SIT8194.0 TRIETHOXYSILYLUNDECANAL C ₁₇ H ₃₆ O ₄ Si coupling agent for DNA HMIS: 2-2-1-X	332.56	150-5°/0.5		1.4343
	5.0g/\$110.00				
	SIT8194.5 TRIETHOXYSILYLUNDECANAL, ETHYLENE GLYCOL ACETAL C ₁₉ H ₄₀ O ₅ Si HMIS: 2-2-1-X	366.60	160-5°/0.25		
	5.0g/\$110.00				
	SIT8185.3 TRIETHOXSILYLBUTYRALDEHYDE, tech-90 C ₁₀ H ₂₂ O ₄ Si contains 3-TRIETHOXSILYL-2-METHYLPROPANAL isomer and cyclic siloxy acetal, 2,2,6-TRIETHOXY-1-OXA-2-SILACYCLOHEXANE	234.37	85-7°/1	0.96	1.414
	[88276-92-0] HMIS: 3-3-1-X	10g/\$140.00			

Developmental










Developmental



A variety of composite materials utilizing methacrylate and aminosilanes are used in laser-printers.

Amino Functional Silanes

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
	Monoamine Functional Silanes - Trialkoxy SIA0610.0 3-AMINOPROPYLTRIETHOXY-SILANE C ₉ H ₂₃ NO ₃ Si AMEO, GAPS flashpoint: 104°C (220°F) ΔHvap: 11.8 kcal/mole viscosity: 1.6 cSt. versatile coupling agent effects immobilization of enzymes ¹ . 1. Enzymes, 84, 55915, 1976	221.37	122-3°/30	0.951	1.4225
	[919-30-2] TSCA HMIS: 3-1-1-X 25g/\$10.00 2.0kg/\$96.00 16kg/\$488.00	TOXICITY- oral rat, LD50: 1780mg/kg primary irritation index: 6.50 γc of treated surface: 37.5 dynes/cm specific wetting surface: 353m ² /g vapor pressure, 100°: 10mm	Commercial		
	SIA0611.0 3-AMINOPROPYLTRIMETHOXY-SILANE C ₆ H ₁₇ NO ₃ Si hydrolysis rate vs AMEO (SIA0610.0): 6:1 [13822-56-5] TSCA HMIS: 3-2-1-X	179.29			
	SIA0587.0 4-AMINO-BUTYLTRIETHOXY-SILANE, 95% C ₁₀ H ₂₅ NO ₃ Si [3069-30-5] HMIS: 2-2-1-X	235.40	114-6°/14	0.941 ²⁵	1.4270 ²⁵
	flashpoint: 109°C (225°F) TOXICITY- oral rat, LD50: 1620mg/kg 10g/\$37.00 50g/\$148.00	Developmental			
	SIA0599.0 m-AMINOPHENYLTRIMETHOXY-SILANE, 90% C ₉ H ₁₅ NO ₃ Si contains other isomers [70411-42-6] HMIS: 3-1-1-X				
	SIA0599.1 p-AMINOPHENYLTRIMETHOXY-SILANE, 90% C ₉ H ₁₅ NO ₃ Si contains other isomers coupler for silica-poly(phenyleneterephthalamide) composite films! 1. J. Mark et al, J. Mater. Chem. 7, 259, 1997 [33976-43-1] HMIS: 3-1-1-X	213.31	110-4°/0.6 (60-2°) mp flashpoint: 180°C (356°F)		
	5.0g/\$76.00	Developmental			
	SIA0599.2 AMINOPHENYLTRIMETHOXY-SILANE, mixed isomers typically 60-70% para, 30-40% meta C ₉ H ₁₅ NO ₃ Si for pure isomers, see SIA0559.0, SIA0559.1 [33976-43-1] HMIS: 3-1-1-X				
	SIA0614.0 3-AMINOPROPYLTRIS(METHOXYETHOXY-ETHOXY)-SILANE, 95% C ₁₈ H ₄₁ NO ₉ Si for melt compounding of polyamide composites [87794-64-7] HMIS: 3-2-1-X	443.61	flashpoint: 68°C (155°F)	1.066	1.448
	25g/\$40.00	Developmental			
	SIA0630.0 11-AMINO-UNDECYLTRIETHOXY-SILANE C ₁₇ H ₃₉ NO ₃ Si contains ~5% isomers [116821-45-5] HMIS: 2-2-1-X				
	SIP6928.0 2-(4-PYRIDYLETHYL)TRIETHOXY-SILANE C ₁₃ H ₂₃ NO ₃ Si see also SIT8396.0, SIP6926.4 HMIS: 3-2-1-X	269.43	105°/0.9 amber liquid	1.00	1.4624 ²⁴
	10g/\$112.00	Developmental			

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰		
	SIT8396.0 2-(TRIMETHOXSILYLETHYL)PYRIDINE C ₁₀ H ₁₇ NO ₃ Si [27326-65-4] HMIS: 3-2-1-X	227.33	105°/0.3 flashpoint: >110°C (>230°F) 10g/\$39.00	1.06	1.4755		
	SIT8410.0 N-(3-TRIMETHOXSILYL)PYRROLE C ₁₀ H ₁₉ NO ₃ Si for electrode modification, polypyrrole adhesion! 1. R. Simon et al. J. Am. Chem. Soc. 104, 2031, 1982. [80906-67-8] HMIS: 3-1-1-X	229.35	105-7°/1 flashpoint: >110°C (>230°F) 5.0g/\$86.00	1.017	1.463	Developmental	
	SIA0598.0 3-(m-AMINOPHENOXY)PROPYLTRIMETHOXY-SILANE, 95% C ₁₂ H ₂₁ NO ₄ Si [55648-29-8] HMIS: 3-1-1-X	271.39	amber liquid 10g/\$53.00	1.02	1.495		
Monoamine Functional Silanes - Water-borne							
	SIA0608.0 AMINOPROPYLSILANETRIOL, 22-25% in water C ₃ H ₁₁ NO ₃ Si mainly oligomers internal hydrogen bonding stabilizes solution [29159-37-3] TSCA HMIS: 2-0-0-X	137.21	flashpoint: >110°C (230°F) pH: 10.0-10.5 25g/\$10.00 2.0kg/\$120.00	1.06		18kg/\$495.00	Commercial
Monoamine Functional Silanes - Dialkoxy							
	SIA0605.0 3-AMINOPROPYLMETHYLDIETHOXSILANE C ₈ H ₂₁ NO ₂ Si coupling agent for foundry resins [3179-76-8] TSCA HMIS: 3-2-1-X	191.34	85-8°/8 TOXICITY- oral rat, LD50: 4760mg/kg flashpoint: 85°C (185°F) 25g/\$10.00	0.916	1.4272	2.0kg/\$172.00	
Monoamine Functional Silanes - Monoalkoxy							
	SIA0602.0 3-AMINOPROPYLDIISOPROPYLETHOXY-SILANE C ₁₁ H ₂₇ NOSi forms hydrolytically stable monolayers [17559-36-1] HMIS: 3-2-0-X	217.43	78-80°/0.4 5.0g/\$45.00	0.872 ²⁵	1.4489	25g/\$180.00	Developmental
	SIA0603.0 3-AMINOPROPYLDIMETHYLETHOXSILANE C ₇ H ₁₉ NOSi Δ Hform: 147.6 kcal/mole [18306-79-1] TSCA HMIS: 3-2-1-X	161.32	78-9°/24 flashpoint: 73°C (163°F) 5.0g/\$48.00	0.857 ²⁵	1.427 ²⁵	25g/\$192.00	
Diamine Functional Silanes - Trialkoxy							
	SIA0591.0 N-(2-AMINOETHYL)-3-AMINOPROPYLTRI-METHOXSILANE N-[3-(TRIMETHOXSILYL)PROPYL]ETHYLENEDIAMINE DAMO C ₈ H ₂₂ N ₂ O ₃ Si visc: 6.5 cSt Ce: 0.8 γ _C , treated surface: 36.5 dynes/cm coupling agent for polyamides and polyesters with good film forming properties coupling agent for brass and copper substrates [1760-24-3] TSCA HMIS: 3-1-1-X	226.36	140°/15 TOXICITY- oral rat, LD50: 7460mg/kg flashpoint: 150°C (302°F) specific wetting surface: 358 m ² /g 25g/\$10.00 2.0kg/\$96.00	1.019 ²⁵	1.450 ²⁵	16kg/\$576.00	Commercial
	SIA0590.5 N-(2-AMINOETHYL)-3-AMINOPROPYLTRI-ETHOXSILANE, 95% C ₁₁ H ₂₈ N ₂ O ₃ Si [5089-72-5] TSCA HMIS: 3-1-1-X	264.5	156°/15 flashpoint: 148°C (298°F) 25g/\$60.00	0.994	1.4367 ²⁵		

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰	
$\text{H}_2\text{NCH}_2\text{CH}_2\text{CH}_2\text{CH}_2\text{CH}_2\text{CH}_2\text{NHCH}_2\text{Si}(\text{OC}_2\text{H}_5)_3$	SIA0592.6 N-(6-AMINOHEXYL)AMINOMETHYL- TRIETHOXYSILANE, 95% $\text{C}_{13}\text{H}_{32}\text{N}_2\text{O}_3\text{Si}$ [15129-36-9] HMIS: 3-2-1-X	292.49	160°/0.1 flashpoint: >110°C (>230°F)	0.928 ²⁵	1.4385 ²⁵	
	25g/\$29.00		100g/\$94.00			
$\text{H}_2\text{N}(\text{CH}_2)_6\text{NHCH}_2\text{CH}_2\text{CH}_2\text{Si}(\text{OCH}_3)_3$	SIA0594.0 N-(6-AMINOHEXYL)AMINOPROPYL- TRIMETHOXYSILANE, 95% $\text{C}_{12}\text{H}_{30}\text{N}_2\text{O}_3\text{Si}$ employed in immobilization of DNA ¹ . immobilizes PCR primers on glass beads ² . 1. C. Kneuer et al, Int'l J. Pharmaceutics, 196(2), 257, 2000. 2. J. Andreadis et al, Nuc. Acid Res., 28, E-5, 2000. [51895-58-0] HMIS: 3-2-1-X	278.47	160-5°/0.15 flashpoint: >110°C (>230°F)	1.11	1.4501	
	10g/\$31.00		50g/\$124.00			
$\text{H}_2\text{NCH}_2\text{CH}_2\text{NH}(\text{CH}_2)_{11}\text{Si}(\text{OCH}_3)_3$	SIA0595.0 N-(2-AMINOETHYL)-11-AMINOUNDECYL- TRIMETHOXYSILANE $\text{C}_{16}\text{H}_{38}\text{N}_2\text{O}_3\text{Si}$ coupling agent with extended spacer-group for remote substrate binding HMIS: 3-1-1-X	334.57	155-9°/0.4 5.0g/\$130.00	0.873 ²⁵	1.4515	
$\text{H}_2\text{NCH}_2\text{CH}_2\text{NHCH}_2$ 	SIA0588.0 (AMINOETHYLAMINOMETHYL)PHENETHYL- TRIMETHOXYSILANE, 90% mixed m,p isomers $\text{C}_{14}\text{H}_{26}\text{N}_2\text{O}_3\text{Si}$ coupling agent for polyimides photochemically sensitive (194nm) ¹ self-assembled monolayers ² . 1. W. Dressick et al, Thin Solid Films, 284, 568, 1996 2. C. Harnett et al, Appl. Phys. Lett., 76, 2466, 2000. [74113-77-2] TSCA HMIS: 3-1-1-X	298.46	126-30°/0.2 flashpoint: > 110°C (>230°F)	1.02	1.5083	
	25g/\$82.00		100g/\$266.00			
$\text{H}_2\text{N}(\text{CH}_2\text{CHO})\text{CH}_2\text{CH}_2\text{NHCH}_2\text{CH}_2\text{CH}_2\text{Si}(\text{OCH}_3)_3$	SIA0599.4 N-3-[(AMINO(POLYPROPYLENOXY))AMINO- PROPYLTRIMETHOXYSILANE 60-65% contains 30-35% amine terminated polypropylene oxide coupling agent with film-forming capability HMIS: 2-2-1-X	337-435	3-4 propyleneoxy units 25g/\$72.00	0.984	1.4508	
Diamine Functional Silanes - Water-borne						
	SIA0590.0 N-(2-AMINOETHYL)-3-AMINOPROPYL- SILANETRIOL, 25% in water mainly oligomers $\text{C}_5\text{H}_{17}\text{N}_2\text{O}_3\text{Si}$ internal hydrogen bonding stabilizes solution [68400-09-9] TSCA HMIS: 2-0-0-X	180.28	flashpoint: >110°C (230°F) pH: 10.0-10.5	1.00		
	100g/\$10.00		2.0kg/\$130.00			
Diamine Functional Silanes - Dialkoxy						
$\text{H}_2\text{NCH}_2\text{CH}_2\text{NHCH}_2\text{CH}_2\text{CH}_2\text{Si}(\text{OCH}_3)_2$	SIA0589.0 N-(2-AMINOETHYL)-3-AMINOPROPYLMETHYL- DIMETHOXYSILANE $\text{C}_8\text{H}_{22}\text{N}_2\text{O}_2\text{Si}$ comonomer for silicones in textile softeners and haircare formulations [3069-29-2] TSCA HMIS: 3-1-1-X	206.36	265° flashpoint: 90°C (194°F) autoignition temp: 280°C specific wetting surface: 380 m ² /g 25g/\$10.00 2.0kg/\$154.00 16kg/\$954.00	0.975 ²⁵	1.4447 ²⁵	
$\text{H}_2\text{NCH}_2\text{CH}_2\text{NHCH}_2\text{CH}(\text{CH}_3)\text{Si}(\text{OCH}_3)_2$	SIA0587.5 N-(2-AMINOETHYL)-3-AMINOISOBUTYL- METHYLDIMETHOXYSILANE, 95% $\text{C}_9\text{H}_{24}\text{N}_2\text{O}_2\text{Si}$ [23410-40-4] TSCA HMIS: 3-2-1-X	220.39	131°/15 flashpoint: 96°C (205°F) 25g/\$90.00	0.960	1.4518	
Diamine Functional Silanes - Monoalkoxy						
$\text{H}_2\text{NCH}_2\text{CH}_2\text{NHCH}_2\text{CH}(\text{CH}_3)\text{Si}(\text{OCH}_3)$	SIA0587.2 (AMINOETHYLAMINO)-3-ISOBUTYLDI- METHYLDIMETHOXYSILANE, 95% $\text{C}_9\text{H}_{24}\text{N}_2\text{OSi}$ [31024-49-4] HMIS: 3-2-1-X	204.39	85-9°/2 25g/\$84.00	0.900 ²⁵	1.4513 ²⁵	

Developmental

Commercial

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
Triamine Functional					
	SIT8398.0 (3-TRIMETHOXYSILYLPROPYL)DIETHYLENE-TRIAMINE, 95% C ₁₀ H ₂₇ N ₃ O ₃ Si hardener, coupling agent for epoxies [35141-30-1] TSCA HMIS: 3-1-1-X	265.43	114-8°/2 flashpoint: 137°C (279°F) γc of treated surface: 37.5 dynes/cm	1.030	1.4590
Secondary Amine Functional					
	SIB1932.2 n-BUTYLAMINOPROPYLTRIMETHOXY-SILANE C ₁₀ H ₂₅ NO ₃ Si coupling agent for urethane coatings [31024-56-3] TSCA HMIS: 2-2-1-X	235.40	102°/3.5 flashpoint: 110°C (230°F)	0.947	1.4246 ²⁵
	SIE4886.0 N-ETHYLAMINOISOBUTYLTRIMETHOXY-SILANE C ₉ H ₂₃ NO ₃ Si adhesion promoter for polyurethane coatings [227085-51-0] TSCA HMIS: 3-2-1-X	221.37	95°/10 flashpoint: 91°C (196°F)	0.952 ²⁵	1.4234
	SIM6500.0 N-METHYLAMINOPROPYLTRIMETHOXY-SILANE C ₇ H ₁₉ NO ₃ Si pK _b ²⁵ H ₂ O: 5.18 orients liquid crystals [3069-25-8] TSCA HMIS: 3-2-1-X	193.32	106°/30 flashpoint: 82°C (179°F) γc of treated surface: 31 dynes/cm	0.978 ²⁵	1.4194
	SIP6724.0 N-PHENYLAMINOPROPYLTRIMETHOXY-SILANE, 95% C ₁₂ H ₂₁ NO ₃ Si oxidatively stable coupling agent for polyimides, phenolics, epoxies [3068-76-6] TSCA HMIS: 3-1-1-X	255.38	132-5°/0.3 flashpoint: 165°C (329°F) specific wetting surface: 307m ² /g	1.07	1.504
	SIA0400.0 3-(N-ALLYLAMINO)PROPYLTRIMETHOXY-SILANE, 95% C ₉ H ₂₁ NO ₃ Si coupling agent for polyesters coupling agent for acrylic coatings for glass containers. 1. Y. Hashimoto et al, Eur. Pat. Appl. EP 289,325, 1988 [31024-46-1] HMIS: 3-2-1-X	219.36	106-9°/25 flashpoint: 88°C (190°F)	0.989 ²⁵	1.4990 ²⁵
	SIC2464.2 (CYCLOHEXYLAMINOMETHYL)TRI-ETHOXY-SILANE, 95% C ₁₃ H ₂₉ NO ₃ Si [26495-91-0] HMIS: 2-1-1-X	275.46	236° flashpoint: 119°C (246°F)	0.95	1.4377
	SIC2464.4 N-CYCLOHEXYLAMINOPROPYLTRIMETHOXY-SILANE C ₁₂ H ₂₇ NO ₃ Si [3068-78-8] HMIS: 3-2-1-X	261.43	114°/3 25g/\$45.00	0.99	1.486 ²⁵
	SIE4885.8 N-ETHYLAMINOISOBUTYLMETHYL-DIETHOXY-SILANE C ₁₁ H ₂₇ NO ₂ Si HMIS: 3-2-1-X	233.43	89°/27 25g/\$72.00		
	SIP6723.67 (PHENYLAMINOMETHYL)METHYL-DIMETHOXY-SILANE, 95% C ₁₀ H ₁₇ NO ₂ Si converts isocyanate terminated polymers to moisture-cureable resins [17890-10-7] HMIS: 3-2-1-X	211.34	255° flashpoint: 106°C (223°F) 25g/\$29.00	1.04	1.5147










Commercial

Developmental

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
	SIP6723.7 N-PHENYLAMINOMETHYLTRIETHOXY-SILANE C ₁₃ H ₂₃ NO ₃ Si [3473-76-5] HMIS: 3-2-1-X	269.42	135-7°/4	1.004 ²⁵	1.485 ²⁵
		25g/\$29.00		100g/\$94.00	
	SIM6498.0 N-METHYLAMINOPROPYLMETHYL-DIMETHOXY-SILANE C ₇ H ₁₉ NO ₂ Si [31024-35-8] HMIS: 3-2-1-X	177.32	93°/25 flashpoint: 80°C (176°F)	0.9173 ²⁵	1.4224 ²⁵
		25g/\$61.00		100g/\$198.00	
Tertiary Amine Functional Silanes					
	SIB1140.0 BIS(2-HYDROXYETHYL)-3-AMINOPROPYL-TRIETHOXY-SILANE, 62% in ethanol C ₁₃ H ₃₁ NO ₅ Si contains 2-3% hydroxyethylaminopropyltriethoxysilane urethane polymer coupling agent employed in surface modification for preparation of oligonucleotide arrays ¹ . 1. G. McGall et al, Proc. Nat'l Acad. Sci., 93, 1355, 1996 [7538-44-5] TSCA HMIS: 3-4-0-X	309.48	flashpoint: 24°C (75°F) specific wetting surface: 252m ² /g	0.92	1.409 ²⁵
		25g/\$30.00		100g/\$98.00	
	SID3395.4 DIETHYLAMINOMETHYLTRIETHOXY-SILANE C ₁₁ H ₂₇ NO ₃ Si catalyst for neutral cure 1-part RTV's [15180-47-9] HMIS: 2-2-1-X	249.43	74-6°/3	0.9336 ²⁵	1.4142 ²⁵
		25g/\$49.00			
	SID3396.0 (N,N-DIETHYL-3-AMINOPROPYL)TRI-METHOXY-SILANE C ₁₀ H ₂₅ NO ₃ Si [41051-80-3] TSCA HMIS: 2-2-1-X	235.40	120°/20 flashpoint: 100°C (212°F)	0.934	1.4245
		25g/\$58.00		100g/\$188.00	
	SID3547.0 3-(N,N-DIMETHYLAMINOPROPYL)TRIMETHOXY-SILANE C ₈ H ₂₁ NO ₃ Si derivatized silica catalyzes Michael reactions ¹ . 1. J. Mdoe et al, Synlett., 625, 1998 [2530-86-1] TSCA HMIS: 2-2-1-X	207.34	106°/30 flashpoint: 99°C (210°F)	0.948 ²⁵	1.4150
		10g/\$30.00		50g/\$120.00	
Quaternary Amine Functional Silanes					
	SIS6994.0 3-(N-STYRYLMETHYL-2-AMINOETHYLAMINO)-PROPYLTRIMETHOXY-SILANE HYDROCHLORIDE, 40% in methanol, inhibited with BHT C ₁₇ H ₃₁ ClN ₂ O ₃ Si see also SIS6993.0 [34937-00-3] TSCA HMIS: 3-4-1-X store <5°	374.98	flashpoint: 13°C (55°F) viscosity: 2.3 cSt	0.91	1.395
		25g/\$10.00		2.0kg/\$244.00	
	SIT8415.0 N-TRIMETHOXY-SILYLPROPYL-N,N,N-TRI-METHYLAMMONIUM CHLORIDE (50% in methanol) N,N,N-TRIMETHYL-3-(TRIMETHOXY-SILYL)-1-PROPANAMINIUM CHLORIDE C ₉ H ₂₄ ClNO ₃ Si employed for bonded chromatographic phases anti-static agent used to treat glass substrates employed in electroblotting [35141-36-7] TSCA HMIS: 2-4-1-X	257.83	flashpoint: 16°C (61°F)	0.927	1.3966
		25g/\$17.00		2.0kg/\$380.00	

Developmental












Commercial

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰	
	SIO6620.0 OCTADECYLDIMETHYL(3-TRIMETHOXYSILYL- PROPYL)AMMONIUM CHLORIDE, 60% in methanol C ₂₆ H ₅₈ ClNO ₃ Si contains 3-5% Cl(CH ₂) ₃ Si(OMe) ₃ employed as lubricant/ anti-static surface treatment orients liquid crystals dispersion/coupling agent for high density magnetic recording media ¹ . application as immobilizeable antimicrobial reported ² . 1. H. Vincent in "Chemically Modified Oxide Surfaces," ed.D. Leyden, Gordon & Breach,1990, p.305 2. W. White et al in "Silanes, Surfaces & Interfaces" ed.D. Leyden, Gordon & Breach, 1986, p.107 [27668-52-6] TSCA HMIS: 3-4-0-X	496.29		0.89		Commercial
	25g/\$24.00	2.0kg/\$280.00				
	SIB0957.0 (2-N-BENZYLAMINOETHYL)-3-AMINOPROPYL- TRIMETHOXYSILANE, hydrochloride 50% in methanol C ₁₅ H ₂₈ N ₂ O ₃ Si.HCl amber liquid [623938-90-9] TSCA HMIS: 3-3-1-X	348.25		0.942	1.4104	Developmental
	25g/\$16.00	100g/\$52.00				
	SID3392.0 N,N-DIDECYL-N-METHYL-N-(3-TRIMETHOXYSILYL- PROPYL)AMMONIUM CHLORIDE, 42% in methanol C ₂₇ H ₆₀ ClNO ₃ Si contains 3-5% Cl(CH ₂) ₃ Si(OMe) ₃ [68959-20-6] TSCA HMIS: 3-4-0-X	510.32		0.863	1.4085	Developmental
	25g/\$46.00					
	SIT7090.0 TETRADECYLDIMETHYL(3-TRIMETHOXYSILYL- PROPYL)AMMONIUM CHLORIDE, 50% in methanol C ₂₂ H ₅₀ ClNO ₃ Si contains 3-5% Cl(CH ₂) ₃ Si(OMe) ₃ [41591-87-1] TSCA HMIS: 3-4-0-X	440.18		0.88	1.397	Developmental
	25g/\$48.00					
	SIT8395.0 N-(TRIMETHOXYSILYLETHYL)BENZYL-N,N,N- TRIMETHYLAMMONIUM CHLORIDE, 60% in methanol C ₁₅ H ₂₈ ClNO ₃ Si candidate for exchange resins and extraction phases HMIS: 3-3-1-X	333.93		0.966		Developmental
	25g/\$80.00					
	SIT8405.0 N-(TRIMETHOXYSILYLPROPYL)ISOTHIO- URONIUM CHLORIDE, 50% in water TRIHYDROXYPROPYLCARBAMIDOTHIOIC ACID HYDROCHLORIDE C ₇ H ₁₉ ClN ₂ O ₃ SSi antimicrobial activity reported [84682-36-0] TSCA HMIS: 2-0-0-X	274.84		1.190	1.441	Developmental
	25g/\$42.00					
Dipodal Amine Functional Silanes						
	SIB1824.5 BIS(TRIETHOXYSILYLPROPYL)AMINE, 95% C ₁₈ H ₄₃ NO ₆ Si ₂ [13497-18-2] TSCA HMIS: 3-1-1-X	425.71	160°/0.6	0.97	1.4265	Commercial
	25g/\$16.00	100g/\$52.00				
	SIB1833.0 BIS(TRIMETHOXYSILYLPROPYL)AMINE, 95% C ₁₂ H ₃₁ NO ₆ Si ₂ dipodal coupling agent [82985-35-1] TSCA HMIS: 3-1-1-X	341.56	152°/4	1.040	1.4320	Commercial
	25g/\$12.00	2.0kg/\$290.00	18kg/\$1170.00			
	SIB1834.0 BIS[(3-TRIMETHOXYSILYL)PROPYL]- ETHYLENEDIAMINE, 62% in methanol C ₁₄ H ₃₆ N ₂ O ₆ Si ₂ dipodal coupling agent for polyamides with enhanced hydrolytic stability provides improved solder resistance for printed circuit boards [68845-16-9] TSCA HMIS: 3-4-1-X	384.62		0.89		Commercial
	25g/\$24.00	2.0kg/\$410.00				

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
	SIB1834.1 BIS[(3-TRIMETHOXYSILYL)PROPYL]- ETHYLENEDIAMINE, 95% C ₁₄ H ₃₆ N ₂ O ₆ Si ₂ coupling agent for polyamides with enhanced hydrolytic stability [68845-16-9] TSCA HMIS: 3-2-1-X	384.62 flashpoint: >110°C (>230°F)		1.050	1.443
	SIB1828.0 BIS[(3-(TRIETHOXYSILYL)PROPYL)UREA, 60% in ethanol C ₁₉ H ₄₄ N ₂ O ₇ Si ₂ [69465-84-5] HMIS: 2-4-1-X	468.73 flashpoint: 24°C (75°F)		0.923	
	SIB1835.5 BIS[(TRIMETHOXYSILYL)PROPYL]UREA, 95% C ₁₃ H ₃₂ N ₂ O ₇ Si ₂ <i>amber liquid</i> flashpoint: >110°C (>230°F) viscosity: 200-250 cSt. [18418-53-6] TSCA HMIS: 3-2-1-X	384.58 flashpoint: >110°C (>230°F) viscosity: 200-250 cSt. 25g/\$19.00		100g/\$62.00	
	SIB1620.0 BIS[(METHYLDIETHOXYSILYL)PROPYL]AMINE, 95% C ₁₆ H ₃₉ NO ₄ Si ₂ dipodal coupling agent [31020-47-0] HMIS: 2-1-1-X	365.66 155°/0.6		0.937	1.4385
	SIB1645.0 BIS[(METHYLDIMETHOXYSILYL)PROPYL]- N-METHYLAMINE, 95% C ₁₃ H ₃₃ NO ₄ Si ₂ HMIS: 3-2-1-X	323.58 140°/2 viscosity: 6-7 cSt. 25g/\$48.00		0.951	1.4368
Specialty Amine Functional Silanes					
	SIT8187.5 N-(3-TRIETHOXYSILYL)PROPYLEN- 4,5-DIHYDROIMIDAZOLE 3-(2-IMIDAZOLIN-1-YL)PROPYLTRIETHOXYSILANE C ₁₂ H ₂₆ N ₂ O ₃ Si coupling agent for elevated temperature cure epoxies utilized in HPLC of metal chelates ¹ . forms proton vacancy conducting polymers w/sulfonamides by sol-gel ² . ligand for molecularly imprinting silica w/ chymotrypsin transition state analog ³ . 1. T. Suzuki et al, Chem. Lett, 881, 1994 2. V. De Zea Bermudez et al, Sol-Gel Optics II, SPIE Proc. 1728, 180, 1992 3. M. Markowitz et al, Langmuir, 1989. [58068-97-6] TSCA HMIS: 2-1-1-X	274.43 134°/2 flashpoint: >110°C (>230°F) viscosity: 5 cSt.		1.005	1.452
	SIU9055.0 UREIDOPROPYLTRIETHOXYSILANE, 50% in methanol C ₁₀ H ₂₄ N ₂ O ₄ Si contains ureidopropyltrimethoxysilane and related transesterification products coupling agent for polyamides, area-formaldehyde resins [23779-32-0] TSCA HMIS: 2-3-1-X	264.40 (-97°)mp flashpoint: 14°C (58°F)		0.92	1.386
	SIA0006.0 ACETAMIDOPROPYLTRIMETHOXYSILANE C ₈ H ₁₉ NO ₄ Si [57757-66-1] HMIS: 3-2-1-X	221.33 162-5°/2-3 10g/\$120.00			1.4410
	SIP6926.2 2-(2-PYRIDYLETHYL)THIOPROPYLTRI- METHOXYSILANE C ₁₃ H ₂₃ NO ₃ SSi chelates metal ions [29098-72-4] HMIS: 3-2-1-X	301.48 156-7°/0.25 10g/\$118.00		1.089	1.498

Developmental


Commercial

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
	SIP6926.4 2-(4-PYRIDYLETHYL)THIOPROPYLTRIMETHOXYSILANE C ₁₃ H ₂₃ NO ₃ Si immobilizeable ligand for immunoglobulin IgG separation using hydrophobic charge induction chromatography (HCIC) [198567-47-4] HMIS: 3-2-1-X	301.48	160-2°/0.2	1.09	1.5037
	pKa: 4.8 10g/\$124.00				
	SID4068.0 3-(1,3-DIMETHYLBUTYLIDENE)AMINO-PROPYLTRIETHOXYSILANE C ₁₅ H ₃₃ NO ₃ Si [116229-43-7] TSCA HMIS: 2-2-1-X	303.52	134°/5 flashpoint: 131°C (268°F) blocked amine - moisture deblocked	0.93	1.437 ²⁵
	25g/\$19.00				
	SIT8394.0 N-[5-(TRIMETHOXYSILYL)-2-AZA-1-OXO-PENTYL]CAPROLACTAM, 95% N-TRIMETHOXYSILYLPROPYLCARBAMOYL CAPROLACTAM patterns in vitro growth of neurons ¹ . 1. J. Hickman et al, J. Vac. Sci. Tech., 12, 607, 1994 C ₁₃ H ₂₆ N ₂ O ₅ Si [106996-32-1] HMIS: 3-1-1-X	318.45	(-39°)mp flashpoint: 136°C (276°F)	1.14	1.4739
	25g/\$24.00				
	SIU9058.0 UREIDOPROPYLTRIMETHOXYSILANE C ₇ H ₁₈ N ₂ O ₄ Si [23843-64-3] TSCA HMIS: 2-3-1-X	222.32	217-225° flashpoint: 99°C (210°F)	1.150	1.386 ²⁵
	25g/\$10.00				
	SID4465.0 N,N-DIOCTYL-N'-TRIETHOXYSILYLPROPYL-UREA C ₂₆ H ₅₆ N ₂ O ₄ Si [259727-10-1] HMIS: 2-2-1-X	488.83		0.924 ²⁵	1.4521 ²⁵
	25g/\$82.00				
Cyclic Azasilanes					
	SIA0380.0 N-ALLYL-AZA-2,2-DIMETHOXYSILACYCLOPENTANE C ₈ H ₁₇ NO ₂ Si [618914-49-1] HMIS: 3-3-1-X	187.31	52-4°/3		
	10g/\$110.00				
	SIA0592.0 N-AMINOETHYL-AZA-2,2,4-TRIMETHYLSILACYCLOPENTANE C ₈ H ₂₁ NSi [18246-33-8] HMIS: 3-2-1-X	156.28	54-6°/2	0.905	1.4768
	10g/\$60.00				
	SIA0604.0 N-(3-AMINOPROPYLDIMETHYLSILA)AZA-2,2-DIMETHYL-2-SILACYCLOPENTANE tech-90 C ₁₀ H ₂₆ N ₂ Si ₂ employed in vapor-phase derivatization of porous sol-gel silica ¹ . 1. D. Brandhuber et al, J. Mater. Chem., 2005 [388606-32-4] HMIS: 3-1-1-X	230.50			1.4705
	10g/\$72.00				
	SIB1932.4 N-n-BUTYL-AZA-2,2-DIMETHOXYSILACYCLOPENTANE C ₉ H ₂₁ NO ₂ Si vapor phase deposition coupling agent for nanoparticles ¹ . 1. B. Arkles et al in "Silanes and Other Coupling Agents, Vol. 3," K. Mittal (Ed.) VSP-Brill, 2004, p179. [618914-44-6] HMIS: 3-2-1-X	203.36	69-71°/3 flashpoint: 85°C (185°F)	0.941	1.438
	25g/\$64.00				
	SID3543.0 2,2-DIMETHOXY-1,6-DIAZA-2-SILACYCLO-OCTANE C ₇ H ₁₈ N ₂ O ₂ Si [182008-07-7] HMIS: 3-2-1-X	190.32	71-3°/2.5 (61-2°)mp		
	25g/\$80.00				
	SIM6501.4 N-METHYL-AZA-2,2,4-TRIMETHYLSILACYCLOPENTANE C ₇ H ₁₇ NSi coupling agent for nanoparticles [18387-19-4] TSCA HMIS: 3-3-1-X	143.30	137°	0.813	1.4308
	25g/\$48.00				

Developmental


© 2006 Gelest, Inc.

Anhydride Functional Silanes




name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
SIT8192.6 3-(TRIETHOXSILYL)PROPYLSUCCINIC ANHYDRIDE, 95% <i>3-(TRIETHOXSILYL)PROPYL)DIHYDRO-3,5-FURANDIONE</i> C ₁₃ H ₂₄ O ₆ Si coupling agent for dibasic surfaces acetic acid-catalyzed hydrolysis yields succinic acid derivative. [93642-68-3] HMIS: 2-1-1-X	304.41	135°/0.2 flashpoint: >100°C (>212°F) viscosity: 20 cSt.	1.070	1.4405
		25g/\$45.00		100g/\$136.00

Azide Functional Silanes




SIA0780.0 6-AZIDOSULFONYLHEXYLTRIETHOXY-SILANE, 95% <i>1-TRIETHOXSILYL-6-SULFONAZIDE-n-HEXANE</i> C ₁₂ H ₂₇ N ₃ O ₅ Si inserts nitrenes into aliphatics and aromatics at temperatures >110°C [96550-26-4] HMIS: 3-2-1-X	353.51	flashpoint: 114°C (237°F)	1.147	1.4634
		25g/\$120.00		


Carboxylate, Phosphonate and Sulfonate Functional Silanes




SIC2263.0 CARBOXYETHYLSILANETRIOL, SODIUM SALT, 25% in water C ₃ H ₆ O ₅ Na ₂ Si [18191-40-7] HMIS: 2-0-0-X	196.14	pH:12-12.5	1.17 ²⁵	
		25g/\$42.00		100g/\$136.00




SIT8189.8 TRIETHOXSILYLPROPYLMALEAMIC ACID, tech 90 C ₁₃ H ₂₅ NO ₆ Si viscosity: 600-900 cSt. may be imidized by heating after deposition [33525-68-7] TSCA HMIS: 3-2-1-X	319.43		1.11	1.472
		25g/\$48.00		



SIT8402.0 N-(TRIMETHOXSILYLPROPYL)ETHYLENE-DIAMINE TRIACETIC ACID, TRISODIUM SALT, 45% in water C ₁₄ H ₂₅ N ₂ Na ₃ O ₉ Si essentially silanetriol, contains NaCl chelates metal ions [128850-89-5] TSCA HMIS: 2-0-0-X	462.42		1.26	
		25g/\$40.00		100g/\$130.00



SIT8378.3 3-(TRIHIDROXSILYL)-1-PROPANE-SULFONIC ACID 30-35% in water C ₃ H ₁₀ O ₆ SSi [70942-24-4] TSCA HMIS: 3-0-0-X	202.26	(-62°)mp pH: <1	1.12	
		25g/\$48.00		












SIC2417.0 2-(4-CHLOROSULFONYLPHENYL)ETHYLTRIMETHOXSILANE, 50% in methylene chloride C ₁₁ H ₁₇ ClO ₅ SSi contains free sulfonic acid; amber color treated silica acts as etherification catalyst ¹ . treatment of surface oxidized PMDSO supports electroosmotic flow ² . 1. B. Sow et al, Microporous & Mesoporous Materials, 79, 129, 2005 2. B. Wang et al, Micro Total Analysis Systems 2004 Vol 2., Roy Soc. Chem., 297, p109 [126519-89-9] HMIS: 3-2-1-X	324.85		1.30 ²⁵	
		25g/\$68.00		100g/\$221.00

Developmental

Masked Carboxylates - See Anhydride and Ester Functional Silanes

PLEASE INQUIRE ABOUT BULK QUANTITIES

Epoxy Functional Silanes

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰	
Epoxy Functional Silanes - Trialkoxy						
	SIE4668.0 2-(3,4-EPOXYCYCLOHEXYL)ETHYL-TRIETHOXYSILANE C ₁₄ H ₂₈ O ₄ Si [10217-34-2] TSCA HMIS: 2-1-1-X	288.46	114-7°/0.4 flashpoint: 104°C (220°F)	1.015	1.4455	Commercial
		25g/\$14.00	100g/\$46.00	2.0kg/\$320.00		
	SIE4670.0 2-(3,4-EPOXYCYCLOHEXYL)ETHYL-TRIMETHOXYSILANE C ₁₁ H ₂₂ O ₄ Si viscosity: 5.2 cSt coefficient of thermal expansion: 0.8 x 10 ⁻³ vapor pressure, 152°: 10mm ring epoxide more reactive than glycidoxypropyl systems. UV initiated polymerization of epoxy group with weak acid donors. forms UV-cureable coating resins by controlled hydrolysis ¹ . 1. J. Crivello et al, Chem. Mater. 9, 1554, 1997.	246.38	95-7°/0.25 TOXICITY- oral rat, LD50: 12,300mg/kg flashpoint: 146°C (295°F) γc of treated surface: 39.5 dynes/cm specific wetting surface: 317 m ² /g	1.065	1.449	Commercial
		[3388-04-3] TSCA HMIS: 3-1-1-X	100g/\$24.00	2.0kg/\$270.00	18kg/\$1044.00	
	SIG5840.0 (3-GLYCIDOXYPROPYL)TRIMETHOXYSILANE 3-(2,3-EPOXYPROPOXY)PROPYLTRIMETHOXYSILANE C ₉ H ₂₀ O ₅ Si coupling agent for epoxy composites employed in electronic "chip" encapsulation.	236.34	120°/2 (<-70°)mp TOXICITY- oral rat, LD50: 8,400 mg/kg	1.070	1.4290	Commercial
		[2530-83-8] TSCA HMIS: 3-1-1-X	100g/\$16.00	2.0kg/\$138.00	18kg/\$693.00	
	SIG5839.0 (3-GLYCIDOXYPROPYL)TRIETHOXYSILANE C ₁₂ H ₂₆ O ₅ Si [2602-34-8] TSCA HMIS:3-2-1-X	278.4	124°/3 flashpoint: 144°C (291°F)	1.00	1.425	Commercial
		[2602-34-8] TSCA HMIS:3-2-1-X	25g/\$38.00	100g/\$124.00	2.0kg/\$580.00	
	SIG5840.1 (3-GLYCIDOXYPROPYL)TRIMETHOXYSILANE 99+% 3-(2,3-EPOXYPROPOXY)PROPYLTRIMETHOXYSILANE C ₉ H ₂₀ O ₅ Si [2530-83-8] TSCA HMIS: 3-1-1-X	236.34	120°/2 (<-70°)mp TOXICITY- oral rat, LD50: 8,400 mg/kg 25g/\$180.00 in fluoropolymer bottle	1.070	1.4290	Developmental
		[2530-83-8] TSCA HMIS: 3-1-1-X				
	SIE4675.0 5,6-EPOXYHEXYLTRIETHOXYSILANE C ₁₂ H ₂₆ O ₄ Si [86138-01-4] HMIS: 3-2-1-X	262.42	115-9°/1.5 flashpoint: 99°C (210°F) 10g/\$84.00	0.960 ²⁵	1.4254 ²⁵	Developmental
		[86138-01-4] HMIS: 3-2-1-X				
Epoxy Functional Silanes - Dialkoxy						
	SIG5832.0 (3-GLYCIDOXYPROPYL)METHYLDIETHOXY-SILANE C ₁₁ H ₂₄ O ₄ Si employed in scratch-resistant coatings for eyeglasses.	248.39	122-6°/5 TOXICITY- oral rat, LD50: >2000mg/kg flashpoint: 122°C (252°F) viscosity: 3.0 cSt	0.978 ²⁵	1.431	Commercial
		[2897-60-1] TSCA HMIS: 2-1-1-X	25g/\$38.00	100g/\$124.00	2.0kg/\$580.00	
	SIG5836.0 (3-GLYCIDOXYPROPYL)METHYLDIMETHOXY-SILANE C ₉ H ₂₀ O ₄ Si relative hydrolysis rate vs. SIG5840.0: 7.5:1	220.34	100°/4 flashpoint: 105°C (221°F)	1.02	1.431 ²⁵	Developmental
		[65799-47-5] TSCA-L HMIS: 3-1-1-X	25g/\$50.00		100g/\$162.00	
Epoxy Functional Silanes - Monoalkoxy						
	SIG5825.0 (3-GLYCIDOXYPROPYL)DIMETHYLETHOXY-SILANE C ₁₀ H ₂₂ O ₃ Si [17963-04-1] TSCA HMIS: 3-2-1-X	218.37	100°/3 flashpoint: 87°C (189°F)	0.950	1.4337 ²⁵	Developmental
		[17963-04-1] TSCA HMIS: 3-2-1-X	10g/\$38.00		50g/\$152.00	

Ester Functional Silanes

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
	SIA0050.0 ACETOXYMETHYLTRIETHOXYSILANE C ₉ H ₂₀ O ₅ Si hydrolyzes to form stable silanol solutions in neutral water [5630-83-1] HMIS: 2-2-1-X	236.34	106°/15	1.042 ²⁵	1.4092
	SIA0055.0 ACETOXYMETHYLTRIMETHOXYSILANE, 95% C ₆ H ₁₄ O ₅ Si [65625-39-0] TSCA-L HMIS: 3-3-1-X	194.26	190-1° flashpoint: 56°C (133°F)	1.085	1.4031
	SIA0100.0 ACETOXYPROPYLTRIMETHOXYSILANE C ₈ H ₁₈ O ₅ Si γc of treated surface: 37.5 dynes/cm [59004-18-1] HMIS: 3-1-1-X	222.31	92°/2 flashpoint: 93°C (200°F)	1.062	1.4146
	SIB0959.0 BENZOYLOXYPROPYLTRIMETHOXYSILANE C ₁₃ H ₂₀ O ₅ Si [76241-02-6] TSCA-L HMIS: 3-2-1-X	284.38	145°/0.2 25g/\$64.00	1.104	1.4806
	SIC2067.0 10-(CARBOMETHOXY)DECYLDIMETHYL-METHOXYSILANE C ₁₅ H ₃₂ O ₃ Si HMIS: 2-1-1-X	288.50	130°/0.3 10g/\$48.00	0.903	1.4399
	SIC2072.0 2-(CARBOMETHOXY)ETHYLTRIMETHOXY-SILANE, 95% contains ~ 20% 1-(carbomethoxy)ethyltrimethoxysilane isomer METHYL(3-TRIMETHOXSILYLPROPIONATE) C ₇ H ₁₆ O ₅ Si [76301-00-3] HMIS: 3-3-1-X	208.29	flashpoint: > 43°C (>110°F) 10g/\$88.00		

Developmental

Halogen Functional Silanes

Halogen Functional Silanes - Trialkoxy

	SIC2295.5 ((CHLOROMETHYL)PHENYLETHYL)-TRIMETHOXYSILANE, mixed m,p isomers C ₁₂ H ₁₉ ClO ₃ Si employed as a high temperature coupling agent ¹ . 1. B. Arkles et al, Modern Plastics, 57(11), 64, 1980.	274.82	115°/1.5 flashpoint: 130°C (282°F)	1.09 ²⁵	1.4930 ²⁵
	SIC2296.2 (p-CHLOROMETHYL)PHENYLTRIMETHOXY-SILANE C ₁₀ H ₁₅ ClO ₃ Si coupling agent for polyimides [24413-04-5] TSCA HMIS: 3-2-1-X	246.77	134-43°/10 flashpoint: 183°C (361°F)	1.14	1.4965
	SIC2298.4 CHLOROMETHYLTRIETHOXYSILANE C ₇ H ₁₇ ClO ₃ Si Grignard reacts with chlorosilanes or intramolecularly to form carbosilanes ¹ . 1. D. Brondani et al, Tet. Lett., 34, 2111, 1993	212.75	90-1°/25 TOXICITY- oral rat, LD50: 2400mg/kg flashpoint: 47°C (117°F)	1.048	1.4069
	SIC2407.0 3-CHLOROPROPYLTRIETHOXYSILANE C ₉ H ₂₁ ClO ₃ Si [5089-70-3] TSCA HMIS: 2-2-0-X	240.80	100-2°/10 flashpoint: 74°C (172°F)	1.009	1.420
	SIC2410.0 3-CHLOROPROPYLTRIMETHOXYSILANE C ₆ H ₁₅ ClO ₃ Si vapor pressure, 100°: 40mm viscosity, 20°C: 0.56 cSt γc of treated surface: 40.5 dynes/cm [2530-87-2] TSCA HMIS: 3-2-1-X	198.72	195-6° flashpoint: 78°C (172°F) TOXICITY- oral rat, LD50: 5628mg/kg specific wetting surface: 394m ² /g	1.077 ²⁵	1.4183 ²⁵

Commercial

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
	SIB1886.0 7-BROMOHEPTYLTRIMETHOXYSILANE C ₁₀ H ₂₃ BrO ₃ Si HMIS: 3-2-1-X	299.28			
			10g/\$86.00		
	SIB1906.0 3-BROMOPROPYLTRIMETHOXYSILANE C ₆ H ₁₅ BrO ₃ Si forms self-assembled monolayers which can be modified w/ pyridine ligands ¹ . 1. S. Paulson et al, J. Chem. Soc., Chem. Commun., 1615, 1992.	243.17	130°/45	1.293	1.440
	[51826-90-5] HMIS: 2-2-1-X		10g/\$64.00		
	SIB1909.0 11-BROMOUNDECYLTRIMETHOXYSILANE, 95% C ₁₄ H ₃₁ BrO ₃ Si contains undecyltrimethoxysilane [17947-99-8] HMIS: 2-1-0-X	355.39	158°/0.8	1.119	1.4559
			10g/\$80.00		
	SIC2298.6 CHLOROMETHYLTRIMETHOXYSILANE C ₄ H ₁₁ ClO ₃ Si [5926-26-1] HMIS: 3-4-1-X	170.67	156° flashpoint: 26°C (79°F)	1.125	1.4070
			10g/\$25.00	50g/\$100.00	
	SII6452.0 3-IODOPROPYLTRIMETHOXYSILANE C ₆ H ₁₅ IO ₃ Si couples zeolite monolayers to glass ¹ 1. K. Ha et al, Adv. Mater., 12(15), 1114, 2002.	290.17	79-80°/2 flashpoint: 78°C (172°F)	1.475	1.4714
	[14867-28-8] HMIS: 3-2-1-X		10g/\$22.00	50g/\$88.00	
	SIT8397.0 3-(TRIMETHOXYSILYLPROPYL)-2-BROMO-2-METHYLPROPIONATE C ₁₀ H ₂₁ BrO ₅ Si for surface initiated ATRP polymerization ¹ . 1. M. Mulvihill et al, J. Am. Chem. Soc., 127, 16040, 2005	329.27	90-5°/0.5	1.243 ²⁵	
	[314021-97-1] HMIS: 2-2-1-X		5.0g/\$180.00		
Halogen Functional Silanes - Dialkoxy					
	SIC2292.0 CHLOROMETHYLMETHYLDIETHOXYSILANE C ₆ H ₁₅ ClO ₂ Si vapor pressure, 70°: 20mm	182.72	160-1°	1.000 ²⁵	1.407
	[2212-10-4] TSCA HMIS: 3-3-1-X		100g/\$110.00		
	SIC2295.2 ((CHLOROMETHYL)PHENYLETHYL)-METHYLDIMETHOXYSILANE mixed m,p isomers C ₁₂ H ₁₉ ClO ₂ Si intermediate for silicone analog of Merrifield resins.	258.82	120-5°/0.5		
	HMIS: 2-1-1-X		25g/\$120.00		
	SIC2355.0 3-CHLOROPROPYLMETHYLDIMETHOXY-SILANE C ₆ H ₁₅ ClO ₂ Si [18171-19-2] TSCA HMIS: 3-2-1-X	182.72	70-2°/11 flashpoint: 80°C (176°F)	1.0250	1.4253
			100g/\$15.00	2.0kg/\$184.00	
Halogen Functional Silanes - Monoalkoxy					
	SIC2278.0 3-CHLOROISOBUTYLDIMETHYLMETHOXY-SILANE C ₇ H ₁₇ ClOSi [18244-08-1] TSCA HMIS: 3-3-1-X	180.75	182° flashpoint: 52°C (125°F)	0.950	1.4331 ²⁵
			25g/\$52.00		
	SIC2286.0 CHLOROMETHYLDIMETHYLETHOXYSILANE C ₅ H ₁₃ ClOSi dipole moment: 2.14 debye	152.70	132-3° flashpoint: 26°C (79°F)	0.944 ²⁵	1.412 ²⁵
	[13508-53-7] TSCA HMIS: 3-3-1-X		25g/\$62.00		
	SIC2337.0 3-CHLOROPROPYLDIMETHYLETHOXYSILANE C ₇ H ₁₇ ClOSi [13508-63-9] HMIS: 2-3-1-X	180.75	87°/30 flashpoint: 46°C (115°F)	0.932 ²⁵	1.4270 ²⁵
			25g/\$48.00		
	SIC2338.0 3-CHLOROPROPYLDIMETHYLMETHOXY-SILANE, 95% C ₆ H ₁₅ ClOSi see also SIC2278.0 [18171-14-7] HMIS: 3-3-1-X	166.73	170-1° flashpoint: 39°C (102°F)	0.9413	1.4278
			10g/\$54.00		

Developmental

Commercial

Developmental

Hydroxyl Functional Silanes

name **MW** **bp/mm (mp)** **D₄²⁰** **n_D²⁰**

Hydroxyl Functional Silanes - Trialkoxy

	SIB1140.0 BIS(2-HYDROXYETHYL)-3-AMINOPROPYL-TRIETHOXYSILANE, 62% in ethanol C₁₃H₃₁NO₅Si contains 2-3% hydroxyethylaminopropyltriethoxysilane urethane polymer coupling agent employed in surface modification for preparation of oligonucleotide arrays ¹ . 1. G. McGall et al, Proc. Nat'l Acad. Sci., 93, 1355, 1996 [7538-44-5] TSCA HMIS: 3-4-0-X	309.48 flashpoint: 24°C (75°F) specific wetting surface: 252m ² /g	0.92	1.409 ²⁵
	25g/\$30.00	100g/\$98.00		

	SIH6172.0 N-(HYDROXYETHYL)-N-METHYLAMINO-PROPYLTRIMETHOXYSILANE, 75% in methanol C₉H₂₃NO₄Si HMIS: 3-3-1-X	237.37 flashpoint: 16°C (61°F) 25g/\$52.00	0.99	1.417
	100g/\$169.00			

	SIH6175.0 HYDROXYMETHYLTRIETHOXYSILANE, 50% in ethanol TRIETHOXYSILYLMETHANOL C₇H₁₈O₄Si contains equilibrium condensation oligomers hydrolysis yields analogs of silica- hydroxymethylsilanetriol polymers ¹ . 1. B. Arkles, US Pat. 5,371,262, 1994 [162781-73-9] HMIS: 2-4-0-X	194.31	0.866	
	25g/\$96.00			

	SIS6995.0 11-(SUCCINIMIDYLOXY)UNDECYL-DIMETHYLETHOXYSILANE, 95% C₁₉H₃₅NO₅Si reagent for immobilization of proteins via primary amines HMIS: 3-2-1-X	385.58	195-200°/0.6 (28°)mp	
	1.0g/\$210.00			

	SIT8192.0 N-(TRIETHOXYSILYL)PROPYL-O-POLY-ETHYLENE OXIDE URETHANE, 95% C₁₀H₂₂NO₄SiO(CH₂CH₂O)₄₋₆H contains some bis(urethane) analog HMIS: 2-1-1-X	400-500 viscosity: 75-125 cSt 25g/\$24.00	1.09	100g/\$78.00
--	--	---	------	--------------

	SIT8189.5 N-(3-TRIETHOXYSILYL)PROPYL-4-HYDROXY-BUTYRAMIDE C₁₃H₂₉NO₅Si anchoring reagent for light directed synthesis of DNA on glass ¹ . 1. G. McGall et al, J. Am. Chem. Soc., 119, 5081, 1997 [186543-03-3] HMIS: 2-2-1-X	307.47 10g/\$29.00		50g/\$116.00
--	---	-----------------------	--	--------------








	SIT8189.0 N-(3-TRIETHOXYSILYL)PROPYL)GLUCONAMIDE 50% in ethanol C₁₅H₃₃NO₉Si water soluble, hydrophilic silane [104275-58-3] HMIS: 2-4-1-X	399.51 flashpoint: 8°C (46°F) 25g/\$26.00	0.951	100g/\$84.00
--	--	---	-------	--------------

	SIB1824.4 2,2-BIS(3-TRIETHOXYSILYL)PROPANOXY-METHYL)BUTANOL, 50% in ethanol C₂₄H₅₄O₉Si₂ for solid state synthesis of oligonucleotides HMIS: 2-4-1-X	542.86	0.899	10g/\$136.00
--	--	--------	-------	--------------






Masked Hydroxyl

	SIT8572.8 11-(TRIMETHYLSILOXY)UNDECYLTRIETHOXY-SILANE C₂₀H₄₆O₄Si₂ masked hydroxyl- deprotected after deposition with acidic aqueous ethanol [75389-03-6] HMIS: 2-1-1-X	406.75 5.0g/\$134.00		
--	---	-------------------------	--	--

Isocyanate and Masked Isocyanate Functional Silanes

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰	
Isocyanate Functional Silanes - Trialkoxy						
	SII6455.0 3-ISOCYANATOPROPYLTRIETHOXY- SILANE, 95% C ₁₀ H ₂₁ NO ₄ Si component in hybrid organic/inorganic urethanes ¹ . 1. S. Cuney et al, Better Ceramics Through Chemistry VII (MRS. Symp. Proc.), 435, 143, 1996 [24801-88-5] TSCA HMIS: 3-2-1-X	247.37	130°/20 flashpoint: 80°C (176°F)	0.99	1.419	Commercial
	25g/\$14.00 100g/\$45.00 2.0kg/\$220.00					
	SII6453.8 (ISOCYANATOMETHYL)METHYLDIMETHOXY- SILANE tech-85 C ₅ H ₁₁ NO ₃ Si reacts w/polymeric diamines to form moisture-cureable polymers [406679-89-8] HMIS: 3-4-1-X store <5°C	161.23	157° flashpoint: 66°C (151°F) autoignition temp.: 290°	1.06	1.435	
	25g/\$82.00					
	SII6456.0 3-ISOCYANATOPROPYLTRIMETHOXY- SILANE, 95% C ₇ H ₁₅ NO ₄ Si [15396-00-6] TSCA HMIS: 3-2-1-X	205.29	95-8°/10 TOXICITY- oral rat, LD50: 878mg/kg viscosity: 1.4 cSt.	1.073	1.4219	
	25g/\$29.00 100g/\$94.00					
Masked Isocyanate						
	SIT8717.0 TRIS(3-TRIMETHOXSILYLPROPYL)ISO- CYANURATE, 95% C ₂₁ H ₄₅ N ₃ O ₁₂ Si ₃ coupling agent for polyimides to silicon metal [26115-70-8] TSCA HMIS: 2-1-1-X	615.86	flashpoint: 102°C (216°F) viscosity: 325-350 cSt.	1.170	1.4610	Commercial
	25g/\$12.00 100g/\$39.00 2.0kg/\$380.00					
	SIT8186.5 (3-TRIETHOXSILYLPROPYL)-t-BUTYL CAR- BAMATE C ₁₄ H ₃₁ NO ₅ Si [137376-38-6] HMIS: 2-1-1-X	321.49	flashpoint: >65°C (>150°F)	0.990	1.4334	Developmental
	25g/\$36.00 100g/\$117.00					
	SIT8188.0 TRIETHOXSILYLPROPYLETHYL CARBAMATE C ₁₂ H ₂₇ NO ₅ Si masked isocyanate [17945-05-0] TSCA HMIS: 2-1-1-X	293.44	124-6°/0.5 flashpoint: 95°C (203°F)	1.015	1.4321	
	25g/\$24.00 100g/\$80.00					
	SIT7908.0 3-THIOCYANATOPROPYLTRIETHOXSILANE C ₁₀ H ₂₁ NO ₃ SSi TOXICITY- oral rat, LD50: 1423mg/kg [34708-08-2] TSCA HMIS: 3-2-1-X	263.43	95°/0.1 flashpoint: 112°C (234°F)	1.03	1.4460	
	50g/\$20.00 250g/\$80.00					









Phosphine and Phosphate Functional Silanes

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
	SIB1091.0 BIS(2-DIPHENYLPHOSPHINOETHYL)- METHYLSILYLETHYLTRIETHOXY-SILANE, mixed isomers C ₃₇ H ₅₀ O ₃ P ₂ Si ₂ analogous structures form ruthenium II complexes w/ high selectivity for hydrogenation ¹ . 1. D. Wu et al, Chem. Mater., 17, 3951, 2005 HMIS: 2-2-1	660.92		1.07	1.5746
	SID4557.5 DIPHENYLPHOSPHINOETHYLDIMETHYL- ETHOXYSILANE C ₁₈ H ₂₅ OPSi [359859-29-3] HMIS: 2-2-1-X	316.46	160°/1	1.004	1.5630
	SID4558.0 2-(DIPHENYLPHOSPHINO)ETHYL- TRIETHOXY-SILANE C ₂₀ H ₂₉ O ₃ PSi immobilizing ligand for precious metals adhesion promoter for gold substrates in microelectronic applications ¹ . forms stable bonds to silica and basic alumina suitable for catalyst immobilization ² . 1. J. Helbert, US Pat. 4,497,890, 1985 2. C. H. Merckle et al, Chem. Mater. 13, 3617, 2001. [18586-39-5] TSCA HMIS: 3-1-0-X	376.50	182°/1.3 flashpoint: 134°C (273°F)	1.05	1.5384
	SID3412.0 DIETHYLPHOSPHATOETHYLTRIETHOXY- SILANE, 95% DIETHOXYPHOSPHORYLETHYLTRIETHOXY-SILANE C ₁₂ H ₂₉ O ₆ PSi water-soluble silane; anti-pilling agent for textiles hydrolysis product catalytically hydrates olefins, forming alcohols ¹ . 1. F. Young et al, US Patent 3,816,550, 1974. [757-44-8] TSCA HMIS: 3-2-1-X	328.41	141°/2 flashpoint: 70°C (158°F)	1.031 ²⁵	1.4216
	SIT8378.5 3-TRIHYDROXYSILYLPROPYLEMETHYL- PHOSPHONATE, SODIUM SALT, 42% in water C ₄ H ₁₂ O ₆ NaPSi contains 4-5% methanol, sodium methylphosphonate [84962-98-1] TSCA HMIS: 1-2-0-X	238.18	flashpoint: 79°C (174°F)	1.25	
		100g/\$15.00		500g/\$60.00	

Developmental

Enabling Your Technology

Sulfur Functional Silanes








	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰		
Sulfur Functional Silanes - Trialkoxy							
	SIM6476.0 3-MERCAPTOPROPYLTRIMETHOXY-SILANE C ₆ H ₁₆ O ₃ SSi	196.34	93°/40	1.051 ²⁵	1.4502 ²⁵	Commercial	
	viscosity: 2 cSt γc of treated surface: 41 dynes/cm specific wetting surface: 348 m ² /g coupling agent for EPDM rubbers and polysulfide adhesives for enzyme immobilization ¹ . treatment of mesoporous silica yield highly efficient heavy metal scavenger ² . employed in coupling of fluorescent biological tags to CdS nanocrystals ³ . 1. Tet. Let., 31, 5773, 1990 2. J. Liu et al, Science, 276, 923, 1997 3. M. Bruohez et al, Science, 281, 2013, 1998.						
	[4420-74-0] TSCA HMIS: 3-2-1-X	100g/\$16.00	2.0kg/\$164.00	18kg/\$630.00			
	SIM6476.1 3-MERCAPTOPROPYLTRIMETHOXY-SILANE 99+%	196.34	93°/40	1.051 ²⁵	1.4502 ²⁵		
	C ₆ H ₁₆ O ₃ SSi low fluorescence grade for high-throughput screening						
	[4420-74-0] TSCA HMIS: 3-2-1-X	25g/\$180.00	in fluoropolymer bottle				
	SIM6475.0 3-MERCAPTOPROPYLTRIETHOXY-SILANE, 95%	238.42	210°	0.9325	1.4331		
	C ₉ H ₂₂ O ₃ SSi flashpoint: 88°C (190°F) TOXICITY- oral rat, LD50: > 2000mg/kg						
	[14814-09-6] TSCA HMIS: 2-2-1-X	25g/\$35.00		100g/\$114.00			
	SID3545.0 2,2-DIMETHOXY-1-THIA-2-SILACYCLO-PENTANE C ₅ H ₁₂ O ₂ SSi	164.29	57-8°/7	1.094			
	reagent for modification of silver and gold surfaces; coupling agent for rubber						
	[26903-85-5] HMIS: 3-3-1-X	25g/\$84.00					
	SIM6480.0 11-MERCAPTOUNDECYLTRIMETHOXY-SILANE	308.55	150°/0.5	0.955	1.4523	Developmental	
	C ₁₄ H ₃₂ O ₃ SSi HMIS: 3-2-1-X						
		2.5g/\$162.00					
	SIO6704.0 S-(OCTANOYL)MERCAPTOPROPYL-TRIETHOXY-SILANE tech-95 C ₁₇ H ₃₆ O ₄ SSi	364.62		0.9686	1.4514		
	masked mercaptan - deblocked w/alcohols latent coupling agent for butadiene rubber						
	[220727-26-4] TSCA HMIS: 2-1-1-X	25g/\$19.00		100g/\$62.00			
	SIP6926.2 2-(2-PYRIDYLETHYL)THIOPROPYLTRI-METHOXY-SILANE C ₁₃ H ₂₃ NO ₃ SSi	301.48	156-7°/0.25	1.089	1.498		
	chelates metal ions						
	[29098-72-4] HMIS: 3-2-1-X	10g/\$118.00					
	SIP6926.4 2-(4-PYRIDYLETHYL)THIOPROPYLTRI-METHOXY-SILANE C ₁₃ H ₂₃ NO ₃ SSi	301.48	160-2°/0.2	1.09	1.5037		
	pKa: 4.8 immobilizeable ligand for immunoglobulin IgG separation using hydrophobic charge induction chromatography (HCIC)						
	[198567-47-4] HMIS: 3-2-1-X	10g/\$124.00					

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰		
	SIT7908.0 3-THIOCYANATOPROPYLTRIETHOXY-SILANE C ₁₀ H ₂₁ NO ₃ SSi NCSCH ₂ CH ₂ CH ₂ Si(OC ₂ H ₅) ₃ TOXICITY- oral rat, LD50: 1423mg/kg coupling agent for butyl rubber in mechanical applications complexing agent for Ag, Au, Pd, Pt ¹ . 1. T. Schilling et al, Mikrochemica Acta, 124, 235, 1996.	263.43	95°/0.1 flashpoint: 112°C (234°F)	1.03	1.4460		
	[34708-08-2] TSCA HMIS: 3-2-1-X	50g/\$20.00		250g/\$80.00			
	SIT8411.0 2-(3-TRIMETHOXY-SILYL-PROPYLTHIO)- THIOPHENE C ₁₀ H ₁₈ O ₃ S ₂ Si HMIS: 3-2-1-X	278.46					
			10g/\$110.00			Developmental	
Sulfur Functional Silanes - Dialkoxy							
	SIM6473.0 MERCAPTOMETHYLMETHYLDIETHOXY- SILANE, 95% C ₆ H ₁₆ O ₂ SSi HMIS: 3-2-1-X	180.34	60°/10 flashpoint: 58°C (136°F)	0.975	1.4446		
			10g/\$90.00				
	SIM6474.0 3-MERCAPTOPROPYLMETHYLDIMETHOXY- SILANE C ₆ H ₁₆ O ₂ SSi intermediate for silicones in thiol-ene UV cure systems	180.34	96°/30 flashpoint: 93°C (199°F)	1.00	1.4502		
			100g/\$23.00	500g/\$92.00			
	[31001-77-1] TSCA HMIS: 3-2-1-X						
Sulfur Functional Silanes - Dipodal							
	SIB1825.0 BIS[3-(TRIETHOXY-SILYL)PROPYL]- TETRASULFIDE, tech-95 TESPT C ₁₈ H ₄₂ O ₆ S ₄ Si ₂ contains distribution of S _n species: n = 2-10, average 3.8 viscosity: 11.2 cSt coupling agent for "green" tires adhesion promoter for precious metals dipodal coupling agent/ vulcanizing agent for rubbers	538.94	250°d TOXICITY- oral rat, LD50: 16,400mg/kg flashpoint: 91°C (196°F)	1.095	1.49		
			25g/\$12.00	100g/\$39.00	2.0kg/\$180.00	Commercial	
	[40372-72-3] TSCA HMIS: 2-2-1-X						
	SIB1824.6 BIS[3-(TRIETHOXY-SILYL)PROPYL]- DISULFIDE, 90% C ₁₈ H ₄₂ O ₆ S ₂ Si ₂ contains sulfide and tetrasulfide	474.82	flashpoint: 75°C (167°F)	1.025	1.457		
			25g/\$32.00	100g/\$104.00			
	[56706-10-6] TSCA HMIS: 2-2-1-X						
	SIB1820.5 BIS-[m-(2-TRIETHOXY-SILYLETHYL)TOLYL]- POLYSULFIDE tech-85 C ₃₀ H ₅₀ O ₆ S(2-4)Si ₂ dark viscous liquid	627-691	flashpoint: 55°C (132°F)	1.10	1.533		
			25g/\$32.00				
	[TSCA HMIS: 2-2-1-X]						
	SIB1827.0 BIS[3-(TRIETHOXY-SILYL)PROPYL]THIO- UREA tech -90 C ₁₉ H ₄₄ N ₂ O ₆ SSi ₂ forms films on electrodes for determination of mercury ¹ . 1. Y. Guo et al, J. Pharm. Biol. Anal., 19, 175, 1999	484.40					
			25g/\$134.00				
	[69952-89-2] HMIS: 2-1-1-X						










Vinyl and Olefin Functional Silanes

name MW bp/mm (mp) D₄²⁰ n_D²⁰

Vinyl and Olefin Functional Silanes - Trialkoxy

	<p>SIA0540.0 ALLYLTRIMETHOXSILANE C₆H₁₄O₃Si adhesion promoter for vinyl-addition silicones [2551-83-9] TSCA HMIS: 3-2-1-X</p>	<p>162.26 146-8° flashpoint: 46°C (115°F) 10g/\$24.00 50g/\$96.00</p>	<p>0.963²⁵ 1.4036²⁵ 2.0kg/\$720.00</p>
	<p>SIS6993.0 3-(N-STYRYLMETHYL-2-AMINOETHYLAMINO)- PROPYLTRIMETHOXSILANE, 40% in methanol inhibited with BHT C₁₇H₃₀N₂O₃Si coupling agent for epoxy composites, primer for epoxy coatings [34937-00-3] TSCA HMIS: 3-4-1-X store <5°</p>	<p>338.52 flashpoint: 13°C (55°F) 25g/\$19.00 100g/\$62.00</p>	<p>0.871 1.3900 2.0kg/\$292.00</p>
	<p>SIV9098.0 VINYLTRIACETOXSILANE C₈H₁₂O₆Si coefficient of thermal expansion: 1.6 x 10⁻³ derivatization byproduct is acetic acid [4130-08-9] TSCA HMIS: 3-2-1-X store <5°</p>	<p>232.26 112-3°/1 flashpoint: 88°C (190°F) 100g/\$16.00</p>	<p>1.167 1.423 2.0kg/\$158.00</p>
	<p>SIV9112.0 VINYLTRIETHOXSILANE C₈H₁₈O₃Si vapor pressure, 20°: 5mm ΔHvap: 6.8 kcal/mole specific heat: 0.25 cal/g/° dipole moment: 1.69 copolymerization parameters- e,Q: -0.42, 0.028 relative rate of hydrolysis vs SIV9220.0: 0.05 [78-08-0] TSCA HMIS: 1-3-1-X</p>	<p>190.31 160-1° TOXICITY- oral rat, LD50: 22,500mg/kg flashpoint: 44°C (111°F) autoignition temperature: 268°C (514°F) γc of treated surface: 25 dynes/cm ΔHform: -463.5 kcal/mole 25g/\$10.00 2.0kg/\$110.00 16kg/\$472.00</p>	<p>0.903 1.3960</p>
	<p>SIV9209.0 VINYLTRIISOPROPENOXSILANE C₁₁H₁₈O₃Si employed as a crosslinker and in vapor phase; derivatization byproduct is acetone [15332-99-7] TSCA HMIS: 1-3-1-X</p>	<p>226.35 73-5°/12 25g/\$19.00 100g/\$62.00</p>	<p>0.926 1.4373</p>
	<p>SIV9210.0 VINYLTRIISOPROPOXSILANE C₁₁H₂₄O₃Si copolymerization parameters- e,Q: -0.36, 0.031 relative rate of hydrolysis vs SIV9220.0: 0.0015 [18023-33-1] TSCA HMIS: 1-3-1-X</p>	<p>232.39 179-81° flashpoint: 51°C (124°F) vapor pressure, 60°: 4mm 25g/\$16.00</p>	<p>0.8659 1.3961²⁵ 100g/\$52.00</p>
	<p>SIV9220.0 VINYLTRIMETHOXSILANE C₅H₁₂O₃Si viscosity: 0.6 cSt copolymerization parameters- e,Q: -0.38, 0.031 employed in two-stage¹ and one-stage² graft polymerization/ cross-linking for PE. copolymerizes with ethylene to form moisture cross-linkable polymers³. 1. H. Scott US Pat. 3,646,155, 1972 2. P. Swarbrick et al, US Pat. 4,117,195, 1978 3. T. Isaka et al, U.S. Pat. 4,413,066, 1983 [2768-02-7] TSCA HMIS: 3-4-1-X</p>	<p>148.23 123° TOXICITY- oral rat, LD50: 8,000mg/kg flashpoint: 28°C (82°F) autoignition temp: 235° vapor pressure, 20°: 9mm 25g/\$10.00 2.0kg/\$96.00 16kg/\$352.00</p>	<p>0.970 1.3930</p>

Commercial












	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
	SIV9275.0 VINYLTRIS(2-METHOXYETHOXY)SILANE C ₁₁ H ₂₄ O ₆ Si vapor pressure, 108°: 2mm employed in peroxide graft-moisture crosslinking of polyethylene relative rate of hydrolysis vs SIV9220.0: 0.50 [1067-53-4] TSCA HMIS: 3-2-1-X	280.39	284-6°	1.0336 ²⁵	1.4271 ²⁵
	TOXICITY- oral rat, LD50: 2960mg/kg flashpoint: 115°C (239°F)	50g/\$10.00		2.0kg/\$110.00	
	SIV9280.0 VINYLTRIS(METHYLETHYLKETOXIMINO)- SILANE, tech-95 C ₁₄ H ₂₇ N ₃ O ₃ Si neutral cross-linker/ coupling agent for condensation cure silicones; byproduct is methylethylketoximine [2224-33-1] TSCA HMIS: 3-3-1-X	313.47	113°/0.1 (-22°)mp	0.982 ²⁵	
	50g/\$15.00		2.0kg/\$180.00		
	SIA0482.0 ALLYLOXYUNDECYLTRIMETHOXY-SILANE C ₁₇ H ₃₆ O ₄ Si ω-olefin for functional self-assembled monolayers HMIS: 2-1-0-X	332.56	140°/0.5	0.914	1.4415
			5.0g/\$124.00		
	SIA0525.0 ALLYLTRIETHOXY-SILANE 3-(TRIETHOXY-SILYL)-1-PROPENE C ₉ H ₂₀ O ₃ Si vapor pressure, 100°: 50mm [2550-04-1] TSCA HMIS 2-3-1-X	204.34	176°	0.9030	1.4074
	flashpoint 47°C (117°F) dipole moment: 1.79 debye 10g/\$29.00		50g/\$128.00		
	SIB0988.0 [(BICYCLOHEPTENYL)ETHYL]TRIMETHOXY- SILANE, 95% endo/exo isomers C ₁₂ H ₂₂ O ₃ Si [68323-30-8] HMIS: 2-1-1-X	242.39	65°/10		
	flashpoint: >110°C (>230°F)		25g/\$120.00		
	SIB0992.0 5-(BICYCLOHEPTENYL)TRIETHOXY-SILANE NORBORNENYLTRIETHOXY-SILANE C ₁₃ H ₂₄ O ₃ Si coupling agent for norbornadiene rubbers component in low dielectric constant films undergoes ring-opening metathetic polymerization (ROMP) with RuCl ₂ (P(C ₆ H ₅) ₃) ₃ ¹ . 1. E. Finkelstein, 10th Int'l Organosilicon Symp. Proc. P-120, 1993 [18401-43-9] TSCA HMIS: 2-2-1-X	256.42	106-8°/8	0.960	1.4486
	flashpoint: 98°C (208°F)		10g/\$30.00	50g/\$120.00	
	SIB1928.0 BUTENYLTRIETHOXY-SILANE, 95% C ₁₀ H ₂₂ O ₃ Si mixed isomers (mainly 3-butenyl) [57813-67-9] HMIS: 2-2-1-X	218.37	64°/6 (-80°)mp	0.90	
	TOXICITY - oral rat, LD50: >5000mg/kg flashpoint: 73°C (163°F) 25g/\$88.00				
	SIC2282.0 2-(CHLOROMETHYL)ALLYLTRIMETHOXY- SILANE C ₇ H ₁₅ ClO ₃ Si versatile coupling agent [39197-94-9] HMIS: 3-2-1-X	210.73	128°/70	1.09	
	flashpoint: 89°C (192°F)		2.5g/\$60.00		
	SIC2459.5 [2-(3-CYCLOHEXENYL)ETHYL]TRIETHOXY- SILANE C ₁₄ H ₂₈ O ₃ Si [77756-79-7] HMIS: 2-1-1-X	272.46	flashpoint: 120°C (248°F)	0.948	1.444
	10g/\$25.00		50g/\$100.00		

Commercial










Developmental

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
	SIC2460.0 [2-(3-CYCLOHEXENYL)ETHYL]TRIMETHOXY SILANE C ₁₁ H ₂₂ O ₃ Si orients liquid crystals in display devices ¹ . coupling agent for aramid fiber reinforced epoxy ² . 1. Sharp, CA101,81758g; Jap. Pat. JP 58122517, 1983 2. U. Lechner, CA112, 218118x; Germ. Offen. DE 3820971, 1989	230.38	109°/6 flashpoint: 80°C (176°F)	1.02	1.4476
	[67592-36-3] TSCA HMIS: 3-2-1-X	10g/\$21.00		50g/\$84.00	
	SIC2520.0 (3-CYCLOPENTADIENYLPROPYL)TRI- ETHOXY-SILANE - dimer C ₁₄ H ₂₆ O ₃ Si may be cracked to monomer at ~190° at 100mm employed in silica-supported purification of fullerenes ¹ . 1. B. Nie et al, J. Org. Chem., 61, 1870, 1996	270.44	115°/0.5 flashpoint: 100°C (212°F)	0.99	
	[102056-64-4] HMIS: 2-2-1-X	10g/\$60.00			
	SID4610.3 (DIVINYLMETHYLSILYLETHYL)TRIETHOXY-SILANE C ₁₃ H ₂₈ O ₃ Si ₂ HMIS: 2-1-1-X	288.54	79-81°/0.15	0.895	
	[52034-14-7] HMIS: 2-1-1-X	10g/\$48.00			
	SID4618.0 DOCOSENYLTRIETHOXY-SILANE, 95% C ₂₈ H ₅₈ O ₃ Si contains internal isomers forms self-assembled monolayers that can be modified to hydroxyls ¹ . 1. J. Peansky et al, Langmuir, 11, 953, 1995	470.88	187-195°/0.05		
	[330457-44-8] HMIS: 1-1-0-X	1.0g/\$110.00			
	SIH5919.0 HEXADECAFLUORODODEC-11-ENYL-1- TRIMETHOXY-SILANE C ₁₅ H ₁₆ O ₃ F ₁₆ Si forms self-assembled monolayers/ reagent for immobilization of DNA HMIS: 3-1-1-X	576.35	90°/0.5		
	[52034-14-7] HMIS: 2-1-1-X	10g/\$48.00			
	SIH6164.2 HEXENYLTRIETHOXY-SILANE C ₁₂ H ₂₆ O ₃ Si primarily α-olefin	246.43	97°/1 flashpoint: 86°C (187°F)	0.883	1.4185
	[52034-14-7] HMIS: 2-1-1-X	10g/\$48.00			
	SIO6709.0 7-OCTENYLTRIMETHOXY-SILANE, 95% C ₁₁ H ₂₄ O ₃ Si contains 10-20% internal olefin isomers coupling agent for "in situ" polymerization of acrylamide for capillary electrophoresis ¹ 1. A. Cifuentes et al, J. Chromatog. A, 830(2), 423, 1999	232.39	48-9°/0.1 flashpoint: 95°C (203°F)	0.940	1.4305
	[52217-57-9] TSCA HMIS: 3-1-1-X	5g/\$36.00		25g/\$144.00	
	SIP6902.6 O-(PROPARGYLOXY)-N-(TRIETHOXY- SILYLPROPYL)URETHANE, 90% C ₁₃ H ₂₅ NO ₅ Si HMIS: 2-2-1-X	303.43	110-20°/0.2 flashpoint: 95°C (203°F) inhibited with MEHQ	0.99	1.4461 ²⁵
	[134000-44-5] HMIS: 2-1-1-X	10g/\$60.00		50g/\$240.00	
	SIS6990.0 STYRYLETHYLTRIMETHOXY-SILANE, 95% C ₁₃ H ₂₀ O ₃ Si inhibited with t-butylcatechol mixed m,p isomers and α,β isomers copolymerization parameter, e,Q: -0.880, 1.500	252.38	98°/0.1 flashpoint: 97°C (207°F)	1.02	1.505
	[134000-44-5] HMIS: 2-1-1-X	10g/\$60.00		50g/\$240.00	
	SIU9049.0 10-UNDECENYLTRIMETHOXY-SILANE C ₁₄ H ₃₀ O ₃ Si HMIS: 2-1-1-X	274.48	102-5°/1	0.908	
	[52217-57-9] TSCA HMIS: 3-1-1-X	5g/\$36.00		25g/\$144.00	
	SIV9088.4 O-(VINYLOXYBUTYL)-N-(TRIETHOXY-SILYL- PROPYL)URETHANE, 95% inhibited w/ MEHQ C ₁₆ H ₂₃ NO ₆ Si UV reactive coupling agent HMIS: 3-2-1-X	363.53		1.015	1.4454
	[52217-57-9] TSCA HMIS: 3-1-1-X	10g/\$82.00			

Developmental

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰	
	SIV9099.0 VINYLTRIS- <i>t</i> -BUTOXYSILANE C ₁₄ H ₃₀ O ₃ Si [5356-88-7] HMIS: 2-2-1-X	274.47	54°/2 flashpoint: 79°C (174°F) 10g/\$48.00	0.869		
	SIV9277.0 VINYLTRIS(METHOXYPROPOXY)SILANE C ₁₄ H ₃₀ O ₆ Si [303746-21-6] HMIS: 2-1-1-X	322.47	flashpoint: 122° (252°F) 25g/\$24.00	0.981	1.424	
Vinyl and Olefin Functional - Dialkoxy						
	SIV9085.0 VINYL METHYLDIETHOXY SILANE C ₇ H ₁₆ O ₂ Si copolymerization parameters- e,Q; -0.86, 0.020 dipole moment: 1.27 [5507-44-8] TSCA HMIS: 2-4-1-X	160.29	133-4° flashpoint: 16°C (61°F) 25g/\$14.00	0.858	1.4000	Commercial
	SIV9086.0 VINYL METHYLDIMETHOXY SILANE C ₅ H ₁₂ O ₂ Si viscosity, 20°: 0.51 cSt [16753-62-1] TSCA HMIS: 3-4-1-X	132.23	104° flashpoint: 8°C (46°F) 25g/\$18.00	0.889	1.395	
Vinyl and Olefin Functional - Monoalkoxy						
	SID4612.0 1,3-DIVINYLTETRAMETHYLDISILAZANE C ₈ H ₁₉ NSi ₂ derivatization byproduct is ammonia adhesion promoter for negative photoresists for silylation of glass capillary columns ¹ 1. M. Jaroniec et al, J. High Resol. Chromatog, 5, 3, 1982 [7691-02-3] TSCA HMIS: 3-3-1-X	185.42	160-1° flashpoint: 34°C (93°F) 50g/\$24.00	0.819	1.4405	
	SIV9072.0 VINYL DIMETHYLETHOXY SILANE C ₆ H ₁₄ O ₂ Si dipole moment: 1.23 [5356-83-2] TSCA HMIS: 2-4-1-X	130.26	99-100° flashpoint: 4°C (39°F) 10g/\$22.00	0.790	1.3983	
	SIT8732.0 TRIVINYLMETHOXY SILANE C ₇ H ₁₂ O ₂ Si HMIS: 3-4-1-X	140.25	131-3° 2.5g/\$14.00		1.4400	Developmental
	SIB1818.0 BIS(TRIETHOXY SILETHYL) VINYL METHYL- SILANE C ₁₉ H ₄₄ O ₆ Si ₃ HMIS: 2-1-1-X	452.82	141-3°/0.15 5.0g/\$120.00	0.943		
	SIB1820.0 BIS(TRIETHOXY SILYL) ETHYLENE, 95% 4,4,7,7-TETRAETHOXY-3,8-DIOXA-4,7-DISILADEC-5-ENE C ₁₄ H ₃₂ O ₆ Si ₂ ~80% trans isomer; contains 1,1-isomer [87061-56-1] HMIS: 2-2-1-X	352.57	122-5°/4 5.0g/\$39.00	0.958	1.4168	
	SIB1832.5 BIS(TRIMETHOXY SILYL METHYL) ETHYLENE C ₁₀ H ₂₄ O ₆ Si ₂ [143727-20-2] HMIS: 3-3-1-X	296.47	5.0g/\$170.00			
	SIB1824.9 1,3-[BIS(3-TRIETHOXY SILYL PROPYL) POLY- ETHYLENOXY]-2-METHYLENE PROPANE C ₅₀ H ₁₀₄ O ₂₀ Si ₂ (average) vinyl functional hydrophilic dipodal coupling agent for protein immobilization HMIS: 2-2-1-X	1113.5	1.0g/\$275.00			

Non-Functional Dipodal Silanes




	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
	SIB1817.0 BIS(TRIETHOXSILYL)ETHANE HEXAETHOXYDISILETHYLENE C ₁₄ H ₃₄ O ₆ Si ₂ ΔHvap: 101.5 kJ/mole additive to silane coupling agent formulations that enhances hydrolytic stability employed in corrosion resistant coatings/primers for steel and aluminum ^{1,2} . sol-gels of α,ω-bis(trimethoxysilyl)alkanes reported ³ . forms mesoporous, derivatizable molecular sieves ⁴ . 1. W. Van Ooij et al, J. Adhes. Sci. Tech. 11, 29, 1997 2. W. Van Ooij et al, Chemtech., 28, 26, 1998. 3. D. A. Loy et al, J. Am. Chem. Soc., 121, 5413, 1999. 4. B. Molde et al, Chem. Mat., 11, 3302, 1999.	354.59	96°/0.3	0.957	1.4052
	[16068-37-4] TSCA-S HMIS: 3-1-1-X	25g/\$15.00	100g/\$49.00	2.0kg/\$420.00	
	SIB1829.0 1,2-BIS(TRIMETHOXSILYL)DECANE C ₁₆ H ₃₈ O ₆ Si ₂ pendant dipodal silane HMIS: 3-2-1-X	382.65	130-2°/0.4	0.984	1.4303
	[16068-37-4] TSCA-S HMIS: 3-1-1-X	25g/\$48.00	100g/\$156.00		
	SIB1830.0 BIS(TRIMETHOXSILYL)ETHANE C ₈ H ₂₂ O ₆ Si ₂ CAUTION: INHALATION HAZARD employed in fabrication of multilayer printed circuit boards 1. J. Palladino, U.S. Pat. 5,073,456, 1991.	270.43	103-4°/5	1.068	1.4091
	[18406-41-2] TSCA HMIS: 4-2-1-X	25g/\$69.00	100g/\$224.00		
	SIB1821.0 BIS(TRIETHOXSILYL)METHANE 4,4,6,6-TETRAETHOXY-3,7-DIOXA-4,6-DISILANONANE C ₁₃ H ₃₂ O ₆ Si ₂ intermediate for sol-gel coatings, hybrid inorganic-organic polymers	340.56	114-5°/3.5	0.9741	1.4098
	[18418-72-9] HMIS: 2-3-0-X	5.0g/\$37.00	25g/\$148.00		
	SIB1832.0 BIS(TRIMETHOXSILYL)HEXANE C ₁₂ H ₃₀ O ₆ Si ₂ sol-gels of α,ω-bis(trimethoxysilyl)alkanes reported ¹ . 1. D. A. Loy et al, J. Am. Chem. Soc., 121, 5413, 1999.	326.54	161°/2	1.014	1.4213
	[87135-01-1] HMIS: 3-2-1-X	10g/\$36.00	50g/\$144.00		
	SIB1824.0 BIS(TRIETHOXSILYL)OCTANE C ₂₀ H ₄₆ O ₆ Si ₂ employed in sol-gel synthesis of mesoporous structures	438.76	172-5°/0.75	0.926	1.4240
	[52217-60-4] TSCA-L HMIS: 2-1-1-X	25g/\$30.00	100g/\$98.00		
	SIB1831.0 BIS(TRIMETHOXSILYLETHYL)BENZENE C ₁₆ H ₃₀ O ₆ Si ₂ mixture of m,p isomers flashpoint: 193°C (380°F)	374.58	148-50°/.1	1.08	1.4734
	[58298-01-4] TSCA HMIS: 2-1-0-X	10g/\$33.00	50g/\$132.00		
	SIT8185.8 1-(TRIETHOXSILYL)-2-(DIETHOXYMETHYL)-SILYLETHANE C ₁₃ H ₃₂ O ₅ Si [18418-54-7] TSCA HMIS: 3-2-1-X	324.56	100°/0.5	0.946	1.4112
	[18418-54-7] TSCA HMIS: 3-2-1-X	25g/\$40.00	100g/\$130.00		
	SIB1660.0 BIS[(3-METHYLDIMETHOXSILYL)PROPYL]-POLYPROPYLENE OXIDE viscosity: 6000-10,000 cSt. w/tin catalyst forms moisture-crosslinkable resins hydrophilic dipodal silane	600-800	flashpoint: >110°C (>230°F)	1.00	
	[75009-88-0] TSCA HMIS: 3-1-1-X	100g/\$24.00	2.0kg/\$288.00		

Commercial

UV Active and Fluorescent Silanes

	name	MW	bp/mm (mp)	n_D^{20}	
	SIB1824.8 BIS(4-TRIETHOXYSILYLPROPOXY-3-METHOXY-PHENYL)-1,6-HEPTADIENE-3,5-DIONE tech-90 $C_{39}H_{60}O_{12}Si_2$ UV: 220, 232(max), 354(broad) metal chelating chromophore HMIS: 2-1-1-X	777.07			
	SID4352.0 3-(2,4-DINITROPHENYLAMINO)PROPYL-TRIETHOXYSILANE, 95% <i>N</i> -[3-(TRIETHOXYSILYL)PROPYL]-2,4-DINITROPHENYLAMINE $C_{15}H_{25}N_3O_7Si$ viscous liquid or solid UV: 222, 258, 350(max), 410 forms χ^2 non-linear optical sol-gel materials by corona poling ^{1,2} . 1. E. Toussaere et al, Non-Linear Optics, 2, 37, 1992 2. B. Lebeau et al, J. Mater. Chem., 4, 1855, 1994 [71783-41-0] HMIS: 2-1-0-X	387.46	(27-30°)mp	1.5665	
	SIH6198.0 2-HYDROXY-4-(3-METHYLDIETHOXYSILYL-PROPOXY)DIPHENYLKETONE, 95% $C_{21}H_{28}O_5Si$ monomer for UV opaque fluids HMIS: 2-1-1-X	388.54			
	SIH6200.0 2-HYDROXY-4-(3-TRIETHOXYSILYLPROPOXY)-DIPHENYLKETONE, 95% $C_{22}H_{30}O_6Si$ density: 1.12 strong UV blocking agent for optically clear coatings, absorbs from 210-420nm UV blocking agent ¹ . 1. B. Anthony, US Pat. 4,495,360, 1985 [79876-59-8] TSCA HMIS: 2-1-1-X	418.56		1.545 ²⁵	
	SIM6502.0 O-4-METHYLCOUMARINYL-N-[3-(TRIETHOXY-SILYL)PROPYL]CARBAMATE $C_{20}H_{29}NO_7Si$ immobilizeable fluorescent compound ¹ . 1. B. Arkles, US Pat. 4,918,200, 1990 [129119-78-4] HMIS: 2-2-1-X	423.54	(88-90°)mp		
	SIT8186.2 7-TRIETHOXYSILYLPROPOXY-5-HYDROXY-FLAVONE, 50% in xylene $C_{24}H_{30}O_7Si$ contains non-reactive dyestuffs HMIS: 2-1-1-X	458.58			
	SIT8187.0 N-(TRIETHOXYSILYLPROPYL)DANSYLAMIDE 5-DIMETHYLAMINO-N-(3-TRIETHOXYSILYLPROPYL)-NAPHTHALENE-1-SULFONAMIDE $C_{21}H_{34}N_2O_5Si$ density: 1.12 fluorescent- employed as a tracer in UV cure composites fluorescence probe for crosslinking in silicones ¹ . 1. P. Leezenberg et al, Chem. Mat., 7, 1784, 1995 [70880-05-6] TSCA HMIS: 2-1-1-X	454.66	115-9°/0.1	1.5421	
	SIT8188.8 2-(2-TRIETHOXYSILYLPROPOXY-5-METHYL-PHENYL)BENZOTRIAZOLE $C_{22}H_{31}N_3O_4Si$ UV blocking agent/stabilizer HMIS: 2-1-1-X	429.59			

Biomolecular Probes

	name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
	SIA0120.0 (N-ACETYLGLYCYL)-3-AMINOPROPYL- TRIMETHOXY-SILANE C ₁₀ H ₂₁ N ₂ O ₇ Si amino-acid tipped silane HMIS: 3-2-1-X	309.37			
			5.0g/\$152.00		
	SIT7909.7 3-(N-THYMIDYL)PROPIONOXY-3- TRIMETHOXY-SILANE C ₁₄ H ₂₄ N ₂ O ₇ Si derivatized surfaces bind adenine modified polymers ¹ . 1. K. Viswanathan et al, Polymer Preprints, 46(2), 1133, 2005 HMIS: 2-2-1-X	360.74			
			1.0g/\$210.00		
	SIT8012.0 DL-α-TOCOPHEROLOXY-3-TRIMETHOXY- SILANE tech-90 C ₃₆ H ₆₄ O ₅ Si HMIS: 2-2-1-X	604.99		0.956	1.485
			10g/\$120.00		

Silyl Hydrides


Silyl Hydrides are a distinct class of silanes that behave and react very differently than conventional silane coupling agents. Their application is limited to deposition on metals (see discussion on p. 17). They liberate hydrogen on reaction and should be handled with appropriate caution.

	SID4629.6 DODECYLSILANE C ₁₂ H ₂₈ Si forms SAMS on gold surfaces 872-19-5 HMIS: 2-2-1-X	200.44	80°/7	0.7753	1.4380 ²⁵
					10g/\$78.00
	SIO6635.0 n-OCTADECYLSILANE C ₁₈ H ₄₀ Si contains 4-6% C ₁₈ isomers forms self-assembled monolayers on titanium ¹ . flashpoint: >110°C (>230°F) 1. A. Fadeau et al, J. Am. Chem. Soc., 121, 12184, 1999 [18623-11-5] TSCA HMIS: 2-1-1-X	284.60	195°/15 (29°)mp	0.794	
					25g/\$46.00
					100g/\$150.00
	SIT8173.0 (TRIDecaFLUORO-1,1,2,2-TETRA- HYDROOCTYL)SILANE C ₈ H ₇ F ₁₃ Si provides vapor-phase hydrophobic surfaces on titanium, gold, silicon [469904-32-3] HMIS: 3-3-1-X	378.22	75°/251.446	1.3184	
					10g/\$190.00
	SIU9048.0 10-UNDECENYLSILANE C ₁₁ H ₂₄ Si HMIS: 2-3-1-X	184.40		0.78	
					2.5g/\$180.00


Organosilane-Modified Silica Nanoparticles

A range of silica structures from 20nm to 1 micron have been modified with silanes to reduce hydroxyl content allowing improved dispersion. Other versions have monolayers with isolated secondary amine functionality, providing controlled



interactions with resins. Systems that maintain low levels of hydroxyls have improved electrical properties. Introduction of low levels of secondary amines impart improved mechanical properties particularly in high humidity environments.





name	MW	bp/mm (mp)	D ₄ ²⁰	n _D ²⁰
SIS6960.0 SILICON DIOXIDE, amorphous fumed silica SiO ₂ surface area, 200m ² /g isoelectric point: 2.2 [112945-52-5] TSCA HMIS: 2-0-0-X	60.09	(>1600°)mp	2.2	1.46




SIS6962.0 SILICON DIOXIDE, amorphous HEXAMETHYLDISILAZANE TREATED fumed silica, HMDZ TREATED SiO ₂ carbon content: 3% approximate ratio: (CH ₃) ₃ Si/HO-Si: 2/1 [68909-20-6] TSCA HMIS: 2-0-0-X	60.09	(>1600°)mp	2.2	1.45
--	-------	------------	-----	------



SIS6962.1M30 SILICON DIOXIDE, amorphous HEXAMETHYLDISILAZANE TREATED fumed silica, HMDZ TREATED SiO ₂ carbon content: 2-3% calculated ratio: (CH ₃) ₃ Si/HO-Si: 1/1 [68909-20-6] TSCA HMIS: 2-0-0-X	60.09	(>1600°)mp	2.2	1.45
--	-------	------------	-----	------



SIS6962.1N30 SILICON DIOXIDE, amorphous CYCLIC AZASILANE/HEXAMETHYLDISILAZANE TREATED fumed silica, N-Methylaminopropylfunctional SiO ₂ carbon content: 4-7% calculated ratio: CH ₃ NHCH ₂ CH ₂ CH ₂ Si/(CH ₃) ₃ Si:HO-Si:1/2/1 [68909-20-6] TSCA HMIS: 2-0-0-X	60.09	(>1600°)mp	2.2	1.45
--	-------	------------	-----	------



Commercial

Developmental

Gelest provides custom surface treatment services. We can handle a wide range of materials with special process considerations including: inert atmospheres, highly flammable and corrosive treatments, as well as thermal and vacuum drying.

PLEASE INQUIRE ABOUT BULK QUANTITIES

Surface Modification with Silanes: What's not covered in "Silane Coupling Agents"?

Polar, hydrophilic and water-dispersible silanes, although important in surface modification, do not have organic functionality and are not discussed with coupling agents. The Gelest brochure entitled "**Hydrophobicity, Hydrophilicity and Silane Surface Modification**" includes these materials.

Chlorosilane, silazane and dialkylaminosilane coupling agents are not discussed in this brochure. These materials can be found in the Gelest catalog entitled "Silanes, Silicones and Metal-Organics." The use of these materials is limited commercially due to the difficulty in handling the corrosive, flammable or toxic byproducts associated with hydrolysis.

Alkyl-silanes and **Aryl-silanes** including **Fluorinated Alkyl-silanes** are important in control of hydrophobicity and surface properties. These materials are discussed in the Gelest brochure "Alkyl-silanes and Aryl-silanes."

Further Reading

Silane Coupling Agents - General References and Proceedings

1. B. Arkles, Tailoring Surfaces with Silanes, CHEMTECH, 7, 766-778, 1977
2. E. Plueddemann, "Silane Coupling Agents," Plenum, 1982.
3. K. Mittal, "Silanes and Other Coupling Agents," VSP, 1992
4. D. Leyden and W. Collins, "Silylated Surfaces," Gordon & Breach, 1980.
5. D. E. Leyden, "Silanes, Surfaces and Interfaces," Gordon & Breach 1985.
6. J. Steinmetz and H. Mottola, "Chemically Modified Surfaces," Elsevier, 1992.
7. J. Blitz and C. Little, "Fundamental & Applied Aspects of Chemically Modified Surfaces," Royal Society of Chemistry, 1999.

Substrate Chemistry - General References and Proceedings

8. R. Iler, "The Chemistry of Silica," Wiley, 1979.
9. S. Pantelides, G. Lucovsky, "SiO₂ and Its Interfaces," MRS Proc. 105, 1988.

Product Information		Molecular Weight	Boiling Point/mm (Melting Point)	Refractive Index	Other Physical Properties
Product Code	Product Name			Specific Gravity	
SIA0588.0	(AMINOETHYLAMINOMETHYL)PHENETHYL-298.46 TRIMETHOXYSILANE, 90% mixed m,p isomers		126-30°/0.2 flashpoint: > 110°C (>230°F)	1.02	1.5083
	<chem>C14H26N2O3Si</chem>				
	coupling agent for polyimides photochemically sensitive (194nm) ¹ self-assembled monolayers ² .				
	1. W. Dressick et al, Thin Solid Films, 284, 568, 1996. 2. C Harnett et al, Appl. Phys. Lett., 76, 2466, 2000.				References
	HYDROLYTIC SENSITIVITY: 7 Si-OR reacts slowly with water/moisture				
[74113-77-2]	TSCA HMIS: 3-1-1-X		25g/\$82.00	100g/\$266.00	
CAS#		Hazardous Rating Information (Health-Flammability-Reactivity)			
Indicates Product listed in TSCA Inventory (L = Low Volume Exemption; S = Significant New Use Restriction)					

PRODUCT NAME	PRODUCT#	PAGE	PRODUCT NAME	PRODUCT#	PAGE
ACETOXYMETHYLTRIETHOXSILANE	SIA0050.0	38	3-(2,3-EPOXYPROPOXY)PROPYLTRIMETHOXSILANE	SIG5840.0	37
ACETOXYMETHYLTRIMETHOXSILANE	SIA0055.0	38	5,6-EPOXYHEXYLTRIETHOXSILANE	SI4675.0	37
(3-ACRYLOXYPROPYL)METHYLDIMETHOXSILANE	SIA0198.0	27	N-ETHYLAMINOISOBUTYLTRIMETHOXSILANE	SI4886.0	31
(3-ACRYLOXYPROPYL)TRIMETHOXSILANE	SIA0200.0	26	(3-GLYCIDOXYPROPYL)DIMETHYLETHOXSILANE	SIG5825.0	37
N-(3-ACRYLOXY-2-HYDROXYPROPYL)-3-AMINOPROPYLTRIMETHOXSILANE	SIA0180.0	26	(3-GLYCIDOXYPROPYL)METHYLDIETHOXSILANE	SIG5832.0	37
3-(N-ALLYLAMINO)PROPYLTRIMETHOXSILANE	SIA0400.0	31	(3-GLYCIDOXYPROPYL)METHYLDIMETHOXSILANE	SIG5836.0	37
ALLYLOXYUNDECYLTRIMETHOXSILANE	SIA0482.0	46	(3-GLYCIDOXYPROPYL)TRIETHOXSILANE	SIG5839.0	37
ALLYLTRIETHOXSILANE	SIA0525.0	46	(3-GLYCIDOXYPROPYL)TRIMETHOXSILANE	SIG5840.0	37
ALLYLTRIMETHOXSILANE	SIA0540.0	45	GLYMO	SIG5840.0	37
N-ALLYL-AZA-2,2-DIMETHOXSILACYCLOPENTANE	SIA0380.0	35	HEXAETHOXYDISILETHYLENE	SIB1817.0	50
AMEO	SIA0610.0	28	2-HYDROXY-4-(3-TRIMETHOXSILYLPROPOXY)DIPHENYLKETONE	SIH6200.0	51
4-AMINOBYTTRIETHOXSILANE	SIA0587.0	28	N-(HYDROXYETHYL)-N-METHYLAMINOPROPYLTRIMETHOXSILANE	SIH6172.0	40
(AMINOETHYLAMINO)-3-ISOBUTYLDIMETHYLMETHOXSILANE	SIA0587.2	30	HYDROXYMETHYLTRIMETHOXSILANE	SIH6175.0	40
(AMINOETHYLAMINOMETHYL)PHENETHYLTRIMETHOXSILANE	SIA0588.0	30	3-(2-IMIDAZOLIN-1-YL)PROPYLTRIMETHOXSILANE	SIT8187.5	34
AMINOPHENYLTRIMETHOXSILANE, MIXED ISOMERS	SIA0599.2	28	3-IODOPROPYLTRIMETHOXSILANE	SIH6452.0	39
N-(2-AMINOETHYL)-11-AMINOUNDECYLTRIMETHOXSILANE	SIA0595.0	30	3-ISOCYANATOPROPYLTRIMETHOXSILANE	SIH6455.0	41
N-(2-AMINOETHYL)-3-AMINOISOBUTYLMETHYLDIMETHOXSILANE	SIA0587.5	30	3-ISOCYANTOPROPYLTRIMETHOXSILANE	SIH6456.0	41
N-(2-AMINOETHYL)-3-AMINOPROPYLMETHYLDIMETHOXSILANE	SIA0589.0	30	3-MERCAPTOPROPYLMETHYLDIMETHOXSILANE	SIM6474.0	43
N-(2-AMINOETHYL)-3-AMINOPROPYLSILANETRIOL	SIA0590.0	30	3-MERCAPTOPROPYLTRIMETHOXSILANE	SIM6475.0	43
N-(2-AMINOETHYL)-3-AMINOPROPYLTRIMETHOXSILANE	SIA0591.0	29	3-MERCAPTOPROPYLTRIMETHOXSILANE	SIM6476.0	43
N-(6-AMINOHEXYL)AMINOMETHYLTRIMETHOXSILANE	SIA0526.0	30	O-(METHACRYLOXYETHYL)-N-(TRIMETHOXSILYLPROPYL)URETHANE	SIM6480.8	43
N-(6-AMINOHEXYL)AMINOPROPYLTRIMETHOXSILANE	SIA0594.0	30	MEMO	SIM6487.4	26
m-AMINOPHENYLTRIMETHOXSILANE	SIA0599.0	28	METHACRYLOXYMETHYLTRIMETHOXSILANE	SIM6482.0	26
N-3-[AMINO (POLYPROPYLENOXY)AMINOPROPYLTRIMETHOXSILANE	SIA0599.4	28	METHACRYLOXYMETHYLTRIMETHOXSILANE	SIM6483.0	26
3-AMINOPROPYLDIMETHYLETHOXSILANE	SIA0603.0	28	METHACRYLOXYMETHYLDIETHOXSILANE	SIM6486.4	27
3-(m-AMINOPHENOXY)PROPYLTRIMETHOXSILANE	SIA0598.0	29	METHACRYLOXYPROPYLETHOXSILANE	SIM6486.5	27
N-AMINOETHYL-AZA-2,2-DIMETHYL-4-METHYLSILACYCLOPENTANE	SIA0592.0	35	METHACRYLOXYPROPYLETHYLDIETHOXSILANE	SIM6486.8	27
p-AMINOPHENYLTRIMETHOXSILANE	SIA0599.1	28	METHACRYLOXYPROPYLETHYLDIMETHOXSILANE	SIM6486.9	27
3-AMINOPROPYLDISOPROPYLETHOXSILANE	SIA0602.0	29	METHACRYLOXYPROPYLTRIMETHOXSILANE	SIM6487.3	26
3-AMINOPROPYLDIMETHYLETHOXSILANE	SIA0603.0	29	METHACRYLOXYPROPYLTRIMETHOXSILANE	SIM6487.4	26
3-AMINOPROPYLMETHYLDIETHOXSILANE	SIA0605.0	29	N-(3-METHACRYLOXY-2-HYDROXYPROPYL)-3-AMINOPROPYLTRIMETHOXSILANE	SIM6481.1	26
3-AMINOPROPYLTRIMETHOXSILANE	SIA0610.0	28	N-METHYLAMINOPROPYLMETHYLDIMETHOXSILANE	SIM6498.0	32
3-AMINOPROPYLTRIMETHOXSILANE	SIA0611.0	28	N-METHYLAMINOPROPYLTRIMETHOXSILANE	SIM6500.0	31
3-AMINOPROPYLTRIS(METHOXYETHOXYETHOXY)SILANE	SIA0614.0	28	NORBORNENYLTRIMETHOXSILANE	SIM6992.0	46
AMINOPROPYLSILANETRIOL	SIA0608.0	29	O-4-METHYLCUMARINYLN-3-(TRIMETHOXSILYL)PROPYL]CARBAMATE	SIM6502.0	51
6-AZIDOSULFONYLHEXYLTRIETHOXSILANE	SIA0780.0	36	7-OCTENYLTRIMETHOXSILANE	SI06709.0	47
5-(BICYCLOHEPTENYL)TRIETHOXSILANE	SIB0992.0	46	OCTADECYLDIMETHYL(3-TRIMETHOXSILYLPROPYL)AMMONIUM Cl	SI06620.0	34
1,2-BIS(TRIMETHOXSILYL)DECANE	SIB1829.0	50	N-PHENYLAMINOMETHYLTRIMETHOXSILANE	SIP6723.7	31
BIS(TRIMETHOXSILYLETHYL)BENZENE	SIB1831.0	50	N-PHENYLAMINOPROPYLTRIMETHOXSILANE	SIP6724.0	31
BENZYOXYOPROPYLTRIMETHOXSILANE	SIB0959.0	38	O-(PROPARGYLOXY)-N-(TRIMETHOXSILYLPROPYL)URETHANE	SIP6902.6	47
BIS(2-HYDROXYETHYL)-3-AMINOPROPYLTRIMETHOXSILANE	SIB1140.0	40	2-(4-PYRIDYLETHYL)TRIETHOXSILANE	SIP6928.0	28
BIS(METHYLDIETHOXSILYLPROPYL)AMINE	SIB1620.0	34	STYRYLETHYLTRIMETHOXSILANE	SI06990.0	47
BIS(TRIETHOXSILYL)ETHANE	SIB1817.0	50	3-(N-STYRYLMETHYL-2-AMINOETHYLAMINO)PROPYLTRIMETHOXSILANE	SI06993.0	45
BIS(TRIETHOXSILYL)ETHYLENE	SIB1820.0	48	3-(N-STYRYLMETHYL-2-AMINOETHYLAMINO)PROPYLTRIMETHOXSILANE HCl	SI06994.0	32
BIS(TRIETHOXSILYL)OCTANE	SIB1824.0	50	TETRADECYLDIMETHYL(3-TRIMETHOXSILYLPROPYL)AMMONIUM Cl	SIT7090.0	33
BIS(TRIETHOXSILYL)ETHANE	SIB1830.0	50	4,4,6,6-TETRAETHOXY-3,7-DIOXA-4,6-DISILANONANE	SIB1821.0	50
BIS(TRIETHOXSILYL)HEXANE	SIB1832.0	50	4,4,7,7-TETRAETHOXY-3,8-DIOXA-4,7-DISILADDEC-5-ENE	SIB1820.0	48
BIS(TRIETHOXSILYL)METHANE	SIB1821.0	50	3-THIOCYANATOPROPYLTRIMETHOXSILANE	SIT7908.0	43
BIS(TRIETHOXSILYLPROPYL)AMINE	SIB1833.0	33	TRIETHOXSILYL BUTYRALDEHYDE	SIT8185.3	27
BIS(3-METHYLDIMETHOXSILYL)PROPYLEPOLYPROPYLENE OXIDE	SIB1660.0	50	(3-TRIMETHOXSILYLPROPYL)-t-BUTYL CARBAMATE	SIT8186.5	41
BIS(3-TRIMETHOXSILYL)PROPYLETHYLENEDIAMINE	SIB1834.1	34	1-(TRIMETHOXSILYL)-2-(DIETHOXYMETHYLSILYL)ETHANE	SIT8185.8	50
BIS(3-TRIMETHOXSILYL)PROPYLETHYLENEDIAMINE, 60%	SIB1834.0	33	1-TRIMETHOXSILYL-6-SUFONAZIDE-n-HEXANE	SIA0780.0	36
BIS(3-(TRIMETHOXSILYL)PROPYL)DISULFIDE	SIB1824.6	44	3-(TRIMETHOXSILYL)PROPYLESUCCINIC ANHYDRIDE	SIT8192.6	36
BIS(3-(TRIMETHOXSILYL)PROPYL)TETRASULFIDE	SIB1825.0	44	3-(TRIMETHOXSILYL)PROPYLEDIHYDRO-3,5-FURANDIONE	SIT8192.6	36
BIS(3-(TRIMETHOXSILYL)PROPYL)UREA	SIB1826.0	34	3-(TRIMETHOXSILYL)PROPYL-p-NITROBENZAMIDE	SIT8191.0	51
3-BROMOPROPYLTRIMETHOXSILANE	SIB1906.0	39	7-TRIMETHOXSILYLPROPOXY-5-HYDROXYFLAVONE	SIT8186.2	51
11-BROMOUNDECYLTRIMETHOXSILANE	SIB1909.0	39	N-(3-TRIMETHOXSILYLPROPYL)-4,5-DIHYDROIMIDAZOLE	SIT8187.5	34
BUTENYLTRIMETHOXSILANE	SIB1928.0	46	N-(3-TRIMETHOXSILYLPROPYL)-4-HYDROXYBUTYRAMIDE	SIT8189.0	40
n-BUTYLAMINOPROPYLTRIMETHOXSILANE	SIB1932.2	31	N-(3-TRIMETHOXSILYLPROPYL)GLUCONAMIDE	SIT8189.0	40
N-n-BUTYL-AZA-2,2-DIMETHOXSILACYCLOPENTANE	SIB1932.4	35	N-(TRIMETHOXSILYLPROPYL)DANSYLAMIDE	SIT8187.0	51
10-(CARBOMETHOXY)DECYLDIMETHOXSILANE	SIC2067.0	38	N-(TRIMETHOXSILYLPROPYL)-O-POLYETHYLENE OXIDE URETHANE	SIT8192.0	50
CARBOXYETHYLSILANETRIOL, SODIUM SALT	SIC2263.0	36	N-TRIMETHOXSILYLPROPYL-O-QUININE URETHANE	SIT8192.4	51
((CHLOROMETHYL)PHENYLETHYL)METHYLDIMETHOXSILANE	SIC2295.2	39	TRIMETHOXSILYLMETHANOL	SIH6175.0	40
((CHLOROMETHYL)PHENYLETHYL)TRIMETHOXSILANE	SIC2295.5	38	TRIMETHOXSILYL MODIFIED POLY-1,2-BUTADIENE	SSP-055	49
(p-CHLOROMETHYL)PHENYLTRIMETHOXSILANE	SIC2296.2	38	TRIMETHOXSILYL MODIFIED POLY-1,2-BUTADIENE	SSP-056	49
2-(4-CHLOROSULFONYLPHENYL)ETHYLTRIMETHOXSILANE	SIC2417.0	36	TRIMETHOXSILYLPROPYLETHYL CARBAMATE	SIT8188.0	41
2-(CHLOROMETHYL)ALLYLTRIMETHOXSILANE	SIC2282.0	46	TRIMETHOXSILYLUNDECANAL	SIT8194.0	27
3-CHLOROISOBUTYLDIMETHYLMETHOXSILANE	SIC2278.0	39	TRIMETHOXSILYLUNDECANAL, ETHYLENE GLYCOL ACETAL	SIT8194.5	22
3-CHLOROPROPYLDIMETHYLMETHOXSILANE	SIC2338.0	39	TRIMETHOXSILYL UNDECANAL, ETHYLENE GLYCOL ACETAL	SIT8405.0	33
3-CHLOROPROPYLMETHYLDIMETHOXSILANE	SIC2355.0	39	3-(TRIMETHOXSILYL)-1-PROPANE-SULFONIC ACID	SIT8378.3	36
3-CHLOROPROPYLTRIMETHOXSILANE	SIC2407.0	38	3-TRIMETHOXSILYLPROPYLETHYLPHOSPHONATE, SODIUM SALT	SIT8378.5	36
3-CHLOROPROPYLTRIMETHOXSILANE	SIC2410.0	38	2-(TRIMETHOXSILYLETHYL)PYRIDINE	SIT8396.0	29
CHLOROMETHYLDIMETHYLETHOXSILANE	SIC2286.0	39	N-(3-(TRIMETHOXSILYL)PROPYLETHYLENEDIAMINE	SIA0591.0	29
CHLOROMETHYLMETHYLDIETHOXSILANE	SIC2292.0	39	(3-TRIMETHOXSILYLPROPYL)DIETHYLENETRIAMINE	SIT8398.0	31
CHLOROMETHYLTRIMETHOXSILANE	SIC2298.4	39	N-(3-TRIMETHOXSILYLETHYL)PROPYLEPYRROLE	SIT8410.0	29
CHLOROMETHYLTRIMETHOXSILANE	SIC2298.6	39	N-(TRIMETHOXSILYLETHYL)BENZYL-N,N,N-TRIMETHYLAMMONIUM Cl	SIT8395.0	33
N-CYCLOHEXYLAMINOPROPYLTRIMETHOXSILANE	SIC2464.4	31	N-(TRIMETHOXSILYLPROPYL)ETHYLENEDIAMINE, TRIACETIC ACID, TRISODIUM SALT	SIT8402.0	36
[2-(3-CYCLOHEXYNYL)ETHYL]TRIETHOXSILANE	SIC2459.5	46	N-(TRIMETHOXSILYLPROPYL)ISOTHIOURONIUM CHLORIDE	SIT8405.0	33
[2-(3-CYCLOHEXYNYL)ETHYL]TRIMETHOXSILANE	SIC2460.0	47	N,N,N-TRIMETHYL-3-(TRIMETHOXSILYL)-1-PROPANAMINIUM Cl	SIT8415.0	32
(3-CYCLOPENTADIENYL)PROPYL]TRIETHOXSILANE-dimer	SIC2520.0	47	N-TRIMETHOXSILYLPROPYL-N,N,N-TRIMETHYLAMMONIUM Cl	SIT8415.0	32
DAMO	SIA0591.0	29	TRIMETHOXSILYLPROPYL MODIFIED (POLYETHYLIMINE)	SSP-060	49
N,N-DIDECYL-N-METHYL-N-METHYL-N-(3-TRIMETHOXSILYLPROPYL)AMMONIUM Cl	SID3392.0	33	TRIS(3-TRIMETHOXSILYLPROPYL)ISOCYANURATE	SIT8717.0	41
DIETHOXYMETHYLSILYL MODIFIED POLY-1,2-BUTADIENE	SSP-058	49	TRIVINYLMETHOXSILANE	SIT8732.0	48
DIETHOXYPHOSPHORYLETHYLTRIMETHOXSILANE	SID3412.0	42	10-UNDECENYLSILANE	SIU9048.0	53
DIETHYLAMINOMETHYLTRIMETHOXSILANE	SID3395.4	32	10-UNDECENYLTRIMETHOXSILANE	SIU9049.0	47
(N,N-DIETHYL-3-AMINOPROPYL)TRIMETHOXSILANE	SID3396.0	32	UREIDOPROPYLTRIMETHOXSILANE	SIU9055.0	34
DIETHYLPHOSPHATOETHYLTRIMETHOXSILANE	SID3412.0	42	UREIDOPROPYLTRIMETHOXSILANE	SIU9058.0	34
2,2-DIMETHOXY-1,6-DIAZA-2-SILACYCLOOCTANE	SID3543.0	35	VINYLDIMETHYLETHOXSILANE	SIU9072.0	48
3-(N,N-DIMETHYLAMINOPROPYL)TRIMETHOXSILANE	SID3547.0	32	VINYLMETHYLDIETHOXSILANE	SIU9085.0	48
5-DIMETHYLAMINO-N-(3-TRIMETHOXSILYLPROPYL)NAPHTHALENE-1-SULFONAMIDE	SIT8187.0	51	VINYLMETHYLDIMETHOXSILANE	SIU9086.0	48
DIMETHOXYMETHYLSILYLPROPYL MODIFIED (POLYETHYLIMINE)	SSP-065	49	VINYLTRIACETOXSILANE	SIU9098.0	45
3-(2,4-DINITROPHENYLAMINO)PROPYLTRIMETHOXSILANE	SID4352.0	51	VINYLTRIETHOXSILANE	SIU9112.0	45
N,N-DIOCTYL-N-TRIMETHOXSILYLPROPYLUREA	SID4465.0	35	VINYLTRISOPROPENOXSILANE	SIU9209.0	45
DIPHENYLPHOSPHINOETHYLDIMETHYLETHOXSILANE	SID4557.2	42	VINYLTRISOPROPOXSILANE	SIU9210.0	45
2-(DIPHENYLPHOSPHINO)ETHYLTRIMETHOXSILANE	SID4558.0	42	VINYLTRIMETHOXSILANE	SIU9220.0	45
N-[3-(TRIMETHOXSILYL)PROPYL]-2,4-DINITROPHENYLAMINE	SID4352.0	51	VINYLTRIS(2-METHOXYETHOXY)SILANE	SIU9275.0	46
1,3-DIVINYLTETRAMETHYLDISILAZANE	SID4612.0	48	VINYLTRIS(METHOXYPROPOXY)SILANE	SIU9277.0	48
DOCOSENYLTRIMETHOXSILANE	SID4618.0	47	VINYLTRIS(METHYLETHYLKETOXIMINO)SILANE	SIU9280.0	46
2-(3,4-EPOXYCYCLOHEXYL)ETHYLTRIMETHOXSILANE	SI4668.0	37	VINYLTRI-t-BUTOXSILANE	SIU9099.0	48
2-(3,4-EPOXYCYCLOHEXYL)ETHYLTRIMETHOXSILANE	SI4670.0	37			

Gelest Product Lines



Silicon Compounds: Silanes & Silicones

Detailed chemical properties and reference articles for over 2000 compounds. The 560 page Gelest catalog of silicon and metal-organic chemistry includes scholarly reviews as well as detailed application information. Physical properties, references, structures, CAS numbers as well as HMIS (Hazardous Material Rating Information) of metal-organic and silicon compounds enable chemists to select materials to meet process and performance criteria.



Reactive Silicones - Forging New Polymer Links

The 48 page brochure describes reactive silicones that can be formulated into coatings, membranes, cured rubbers and adhesives for mechanical, optical, electronic and ceramic applications. Information on reactions and cures of silicones as well as physical properties shortens product development time for chemists and engineers. The detailed text provides starting-point formulations, references and application information. Vinyl, hydride, silanol and alkoxy functional silicones are provided for conventional silicone cure systems. Amine, epoxy, methacrylate, hydroxy and mercapto silicones are provided for hybrid organic-silicone cure systems.



Silicone Fluids - Stable, Inert Media

Design and Engineering properties for conventional silicone fluids as well as thermal, fluorosilicone, hydrophilic and low temperature grades are presented in a 24 page selection guide. The brochure provides data on thermal, rheological, electrical, mechanical and optical properties for silicones. Silicone fluids are available in viscosities ranging from 0.65 to 2,500,000 cSt.






Alkyl-Silanes and Aryl-Silanes

A description of non-functional silanes that are used to prepare hydrophobic and water repellent surfaces, specialty resins and modified ceramics is given in an 8 page brochure. The emphasis is on distinguishing the features and benefits of the entire range of commercial alkyl-silanes and aryl-silanes, including fluorinated alkyl-silanes.






Metal-Organics for Material & Polymer Technology




A reference manual for optical and electronic and nanotechnology applications. The literature provides information on metallization, electroceramic, and dielectric applications of silicon, germanium, aluminum, gallium, copper and other metal chemistries. Deposition techniques include ALD, CVD, spin coating and self-assembled monolayers (SAMs). Presents chemistry and physics in the context of device applications ranging from ULSI semiconductors to DNA array devices to flat-panel displays.



Silane
Endcapper



Urethane
Polymer
Backbone



Silane
Endcapper

Gelest

Gelest Inc.
11 East Steel Road
Morrisville, PA 19067
Phone (215) 547-1015
Fax: (215) 547-2484
www.gelest.com